




VIDAS[®]

The pioneer



PIONEERING DIAGNOSTICS

VIDAS®

The pioneer



VIDAS® offers your lab simplicity and productivity, with its loading capacity of 30 tests. It will enhance your workflow, allowing you to handle batch series while maintaining specific management of emergency tests, thanks to its five independent compartments.

A Customized Solution for your Daily Workload

DESIGN


- Streamline your workflow by optimizing test loading
- Increase user-friendliness for data entry and reporting patient results
- PC with external barcode reader & printer

RELIABLE

- System robustness (MTBF⁽¹⁾ > 700 days)
- Single dose concept, ready-to-use reagents
- Guarantees fast, reliable results, processing up to 10 parameters in one run
- Just Load & Go

FLEXIBLE

- Five independent 6-test sections
- Processing of different parameters simultaneously
- Packaging suite to different test volumes (30 or 60 tests)
- 100 parameters available⁽²⁾



	VIDAS®
Weight	65 kg (143 lbs)
Voltage	100-240 VAC
Electric Consumption	3-1.2 A
Frequency	50-60 Hz
Power	280 VA

Discover the VIDAS community and become a member by registering on www.myvidas.com



* SPR: Solid Phase Receptacle
 (1) MTBF = Mean Time Between Failure
 (2) Pending list of assay registered in a specific country