

# We protect you

Medical Products  
Clothing and Gloves  
from own production



**sänger**<sup>®</sup>



# Welcome!

Dear Customer,

***We are happy to present you  
our new Medical Catalogue 2015***

The new catalogue includes a comprehensive range of our medical goods. The descriptions of our products will give you more precise and more detailed information about handling and properties of the goods.

We have improved the functionality of the catalogue to accelerate your search for the most suitable product that corresponds to your requirements.

The Directives and Norms that we are complying to are also listed in the catalogue in order to illustrate our high quality standards.

Our new Webshop is now easier to use. We welcome you to our Website: ***<http://www.sanger.de>***

You can be assured: Sänger GmbH will be a qualified and reliable partner at your side.

We are looking forward to supporting you with your daily work.


Yours sincerely

Sänger GmbH



**sänger**<sup>®</sup>





# Table of Contents

European Directives and Norms	6
National and International Standards	7
Examination Gloves Latex, non-sterile	8-12
Examination Gloves Nitrile, non-sterile	13-15
Examination Gloves Neoprene, non-sterile	15
Examination Gloves Vinyl, non-sterile	16
Examination Gloves Polyethylene, non-sterile	17-18
Examination Gloves Accessories, non-sterile	18-19
Examination Gloves, sterile	20
Surgical Gloves, sterile	21-23
Special Gloves	24-26
Disinfectants	27-29
Skin Care	30-31
Face Masks, Respiratory Protection	32-34
Surgical Caps	35
Clothing, non-sterile	36-40
Easy and Clean sticky mat	40
Clothing, sterile	41
Surgical Retractor Covers	42





Surgical Barriers	43-44
Gynaecological Equipment	45
Softimol Bed Protection, Disposable Bed Sheet	46
Syringes, Cannulas	47-48
Waste Removal Safety Containers	49
Plasters, Sport Tape, Pretape, Cohesive Bandage	50-53
Cooling and Heating Products	54-55
Hot Water Bottles	56
Napkins	56
Medical Accessories	57-60
Wellness	61-62
Anna's Compresses	63
Rubber Ramps	64
Bath and Shower Mats	65
Overview Chemical Resistances	66
Order Form	67

## Directives and Norms

For the general certification and assessment of medical protective equipment in Europe two EU directives are crucial. Following below, you will find the relevant classifications and pictograms. According to each individual certification, you will find the pictograms in the product descriptions.




### 1. European Directive concerning Medical Devices 93/42/EEC

There are a total of six classifications concerning medical products: class I, Is, Im, IIa, IIb and III. If independent testing of the product is required, it is done by a certified institute and the product is marked with the four-digit identification number of the institute. For the products in this catalogue the following classes are relevant:












-  **Class I**  
Low level of invasiveness, temporary use
-  **Class Is**  
Low level of invasiveness, temporary use, sterile; independent evaluation
-  **Class IIa**  
Moderate level of invasiveness, temporary use, sterile; independent evaluation
-  **Class IIb**  
Systemic effects, long-term use; independent evaluation

### 2. European Directive concerning Personal Protective Equipment 89/686/EEC

There are a total of three categories. If independent testing of the product is required, it is done by a certified institute and the product is marked with the four-digit identification number of the institute:

-  **Category I**  
For areas with low risk
-  **Category II**  
For areas with medium risk, EC-model testing
-  **Category III**  
For areas with high, irreversible or fatal risk, EC-model testing

### Additionally protective clothing and gloves are marked with the following pictograms, according to each norm:

- |   |  |
|---|--|
|  General requirements for gloves (according to EN 420)                        |  General Requirements for disposable medical gloves (according to EN 455)               |
|  Instruction manual (according to EN 374)                                     |  Protection against microorganisms (according to EN 374)                                |
|  Protection against chemical hazards (according to EN 374-1:2003, 5.2.1)      |  High protection against chemical hazards (according to EN 374-1:2003, 5.2.1 and 5.3.2) |
|  Protection against ionizing radiation (according to EN 421)                  |  Protection against radioactive particle contamination (according to EN 421)            |
|  Radiation protection in medical diagnostic radiology (according to EN 61331) |  Protection against mechanical hazards (according to EN 388)                            |
|  For use with food  |  |
- |                         |                   |                                 |                   |
|-------------------------|-------------------|---------------------------------|-------------------|
| Criteria                | Performance Level | Criteria                        | Performance Level |
| 1.) Abrasion resistance | 0-4               | 3.) Tear propagation resistance | 0-4               |
| 2.) Cut resistance      | 0-5               | 4.) Puncture resistance         | 0-4               |







## National and International Standards

In order to ensure high quality and satisfied customers, Sanger GmbH complies with the following national and international standards:

- Directive for Medical Devices 93/42/EEC
- Directive for Personal Protective Equipment 89/686/EEC
- DIN EN ISO 9001:2008 Quality Management Systems
- DIN EN ISO 13485:2012 Medical Product Quality Management Systems
- The German Law about Medical Products (MPG)



### Additionally our products are tested according to the following standards:

-  **AQL 1,5**  
AQL (Accepted Quality Level)  
According to Norm EN 455 and ISO 2859-1 there is a quality management requirement on spot tests. It is tested for defects, permeability and physical qualities and there are several quality levels (<0,65; <1,0; <1,5; <2,0; <4,0), with a lower value standing for higher quality.
-  **FDA Standard**  
Tested with the FDA (Food and Drug Administration) 1000 ml water hold test for medical gloves.
-  **OSHA Standard**  
Complies with the OSHA (Occupational Safety and Health Administration) Standard for Hand Protection.
-  **ICRP & ICRU Tested**  
Conform to the suggestions in ICRP (International System of Radiological Protection) Report #60 and ICRU (International Commission on Radiation Units and Measurements) Report #51 concerning radiation protection gloves.
-  **EN 14683**  
Conform to EN 14683 Type II norm for surgical masks.
-  **EN 149:2001**  
Conform to EN 149:2001 Type II norm for fine particle filters.



# Latex

Latex is the milk-like sap from the rubber tree (*Hevea Brasiliensis*). It is tapped from the tree and is refined into a pure latex liquid form via a complex centrifugal process. Because of its polymere structure and its high elasticity, this is the best material for manufacturing medical gloves. Latex gloves are comfortable, flexible, tear-resistant and guarantee a good grip. They provide the best protection against microorganisms like viruses, bacteria etc. All Sanger disposable examination gloves made of latex are a natural product, ambidextrous, thiuram-free, have a low protein content <10 µg/g (Modified Lowry Test, TÜV Rheinland) and a beaded cuff edge.


## Plus for better protection



**PLUS textured**

Size	Art.-No.
extra small (5 - 5½)	01 145
small (6 - 6½)	01 146
medium (7 - 7½)	01 147
large (8 - 8½)	01 148
extra large (9 - 9½)	01 149

**Packing** 100 pieces per dispenser box  
10 dispensers per outer carton



light yellow  
24 cm

The Plus glove is an excellent choice when it comes to medical gloves. Its surface provides a very good grip and good tactile properties for use with instruments or dealing with patients. The glove gives the ideal fit and easy donning.


- Powder Free
- Textured Surface
- Inner and Outer Chlorinated



**PLUS28 textured**

Size	Art.-No.
small (6 - 6½)	01 546
medium (7 - 7½)	01 547
large (8 - 8½)	01 548
extra large (9 - 9½)	01 549

**Packing** 100 pieces per dispenser box  
10 dispensers per outer carton



light yellow  
28 cm

Standards:



## Saturn for use with instruments




A special treatment to the surface of the Saturn glove provides smooth donning and very good grip on instruments, even in a wet environment.

- Powder Free
- Textured Surface
- Inner Chlorinated and Outer Polymer Coated

**SATURN textured**

Size	Art.-No.
extra small (5 - 5½)	01 335
small (6 - 6½)	01 336
medium (7 - 7½)	01 337
large (8 - 8½)	01 338
extra large (9 - 9½)	01 339


**Packing** 100 pieces per dispenser box  
10 dispensers per outer carton



natural  
24 cm



Standards:



## Touch for treating patients


The surface of the Touch glove is specially treated for the contact with patients. It does not pull on hair or skin and the smooth surface guarantees very good tactile properties.



**TOUCH smooth**

Size	Art.-No.
extra small (5 - 5½)	01 820
small (6 - 6½)	01 821
medium (7 - 7½)	01 822
large (8 - 8½)	01 823
extra large (9 - 9½)	01 824


**Packing** 100 pieces per dispenser box  
10 dispensers per outer carton



light yellow  
24 cm

- Powder Free
- Smooth Surface
- Inner and Outer Chlorinated

Standards:





## Profi for optimum grip

- natural rubber latex
- powderfree
- inner and outer polymer coated
- textured surface

Excellent tactile properties with smooth inner surface. Optimum grip make the gloves slip-resistant, even under wet conditions. Gloves offer a secure grip on instruments, glass, etc. Very low protein content of less than 10 µg/g Modified Lowry Test, Dermatological University of Erlangen/ Germany. Hypoallergenic. Chlorine free. For single use.



### PROFI

Size	Art.-No.
extra small (5 - 5½)	01 140
small (6 - 6½)	01 141
medium (7 - 7½)	01 142
large (8 - 8½)	01 143
extra large (9 - 9½)	01 144

Packing: 100 pieces per dispenser box  
10 dispensers per outer carton



#### Standards:



## Chemotex for high risk areas

The extra long and extra thick Chemotex glove is used in areas where maximum protection is required. The structure of the glove is reinforced in order to achieve the highest barrier and tear-resistance values. The textured surface helps to handle objects save and easy.

- Powder Free
- Textured Surface
- Extra Long and Extra Thick
- Inner and Outer Chlorinated




### CHEMOTEX textured

Size	Art.-No.
small (6 - 6½)	01 281
medium (7 - 7½)	01 282
large (8 - 8½)	01 283
extra large (9 - 9½)	01 284

Packing: 50 pieces per dispenser box  
10 dispensers per outer carton



#### Standards:



## Verde for sensitive hands




With soothing Aloe Vera



The Verde glove is a both sides chlorinated, surface-textured glove with an Aloe Vera solution. The Aloe Vera has a healing and skin care effect, being ideal for users with sensitive skin. The glove contains a pleasant Aloe Vera scent.

- Powder Free
- Textured Surface
- Inner Aloe-Vera
- Inner and Outer Chlorinated

#### Standards:



### VERDE textured

Size	Art.-No.
extra small (5 - 5½)	01 645
small (6 - 6½)	01 646
medium (7 - 7½)	01 647
large (8 - 8½)	01 648
extra large (9 - 9½)	01 649

Packing: 100 pieces per dispenser box  
10 dispensers per outer carton



## Sedna against light reflexes



The black coloured Sedna glove reduces light reflections, making it ideal for laser treatments.

- Powder Free
- Textured Surface
- Inner and Outer Chlorinated


### SEDNA textured

Size	Art.-No.
extra small (5 - 5½)	01 745
small (6 - 6½)	01 746
medium (7 - 7½)	01 747
large (8 - 8½)	01 748
extra large (9 - 9½)	01 749

Packing: 100 pieces per dispenser box  
10 dispensers per outer carton



#### Standards:





## Griff *the classic powdered latex glove*



The distinctive element of the economical Griff glove is the textured surface. Its lightly powdered with food-grade corn starch powder which allows easy donning.


- Lightly Powdered
- Textured Surface

### GRIFF *smooth*

Size	Art.-No.
extra small (5 - 5½)	01 130
small (6 - 6½)	01 131
medium (7 - 7½)	01 132
large (8 - 8½)	01 133
extra large (9 - 9½)	01 134
Packing	100 pieces per dispenser box 10 dispensers per outer carton



Standards:



## Finger Cots *for special treatment*

Finger cots are made out of latex, provide a very good fit and reduce hand perspiration. They are an alternative if only finger protection is needed.


- Powder free
- Silicone free
- Latex

### FINGER COTS *powder free*

Size	Art.-No.
small	01 271
medium	01 272
large	01 273
extra large	01 274
Packing:	100 pieces per polybag

### FINGER COTS *powder free*

Size	Art.-No.
small	01 471
medium	01 472
large	01 473
extra large	01 474
Packing:	1440 pieces per polybag



## Nitrile

Nitrile is an elastomer based on petroleum. It is produced through copolymerization, which results in a material that is a good alternative to Latex. Nitrile is the material of choice for users with a Latex allergy and also offers a high resistance against chemicals. All Sängers disposable examination gloves made of nitrile are ambidextrous, thiuram-free and have a beaded cuff edge.

## Alfatex *universally usable*

This Nitrile glove is suited for a variety of situations, and serves to handle chemicals. The textured surface guarantees good grip as well as tactile properties.

- Latex Free
- Powder Free
- Textured Surface




### ALFATEX *textured*

Size	Art.-No.
extra small (5 - 5½)	01 170
small (6 - 6½)	01 171
medium (7 - 7½)	01 172
large (8 - 8½)	01 173
extra large (9 - 9½)	01 174
Packing:	100 pieces per dispenser box 10 dispensers per outer carton



Standards:



### ALFATEX *textured*

Size	Art.-No.
extra small (5 - 5½)	01 370
small (6 - 6½)	01 371
medium (7 - 7½)	01 372
large (8 - 8½)	01 373
extra large (9 - 9½)	01 374
Packing:	100 pieces per dispenser box 10 dispensers per outer carton



## Eris *prevents light reflexes*



### ERIS

Size	Art.-No.
small (6 - 6½)	01 851
medium (7 - 7½)	01 852
large (8 - 8½)	01 853
extra large (9 - 9½)	01 854
Packing:	100 pieces per dispenser box 10 dispensers per outer carton



The black coloured Eris glove reduces light reflections, making it ideal for laser treatments.

- Powder Free
- Textured Surface
- Inner and Outer chlorinated







## Alfatex30 for use with chemicals

This Nitrile glove is universally usable, but is especially suitable for the use with chemicals. Its long cuff ensures better protection against fluids. The textured surface guarantees good grip and tactile properties.

- Latex Free
- Powder Free
- Textured Surface
- Extra long

### ALFATEX 30 textured

Size		Art.-No.
small	(6 - 6½)	01 571
medium	(7 - 7½)	01 572
large	(8 - 8½)	01 573
extra large	(9 - 9½)	01 574

Packing: 100 pieces per dispenser box  
10 dispensers per outer carton



Standards:



## Neoprene


Polychloroprene, better known as Neoprene, is a synthetic rubber made for various purposes. Neoprene gloves have an elasticity and protect against micro-organisms similar to latex but provide a much higher mechanical resistance. Thiuram-free and have a beaded cuff edge.

## Neotril for the best allround protection

The Neotril is produced with the most modern manufacturing processes. A special Neoprene composite material gives this glove the advantages of three materials, while having none of the disadvantages. Neotril possesses the elasticity and protection against biological hazards of latex, the resistance against chemicals of nitrile, and the high tolerance against mechanical stress of neoprene. A textured surface guarantees good grip and good tactile properties. The Neotril is the most powerful existing disposable glove.

- Neoprene Composite Material
- Powder Free
- Inner and Outer Chlorinated
- Textured Surface


Standards:



### NEOTRIL textured

Size		Art.-No.
extra small	(5 - 5½)	01 770
small	(6 - 6½)	01 771
medium	(7 - 7½)	01 772
large	(8 - 8½)	01 773
extra large	(9 - 9½)	01 774

Packing: 100 pieces per dispenser box  
10 dispensers per outer carton



## Betatex for large-scale use

The new Betatex glove provides good protection while being well-priced, ideal for the use in areas where large quantities of gloves are required. The textured surface guarantees good grip and good tactile properties.

- Latex Free
- Powder Free
- Textured Surface



Standards:



### BETATEX textured

Size		Art.-No.
extra small	(5 - 5½)	01 575
small	(6 - 6½)	01 576
medium	(7 - 7½)	01 577
large	(8 - 8½)	01 578
extra large	(9 - 9½)	01 579

Packing: 150 pieces per dispenser box  
10 dispensers per outer carton



### BETATEX textured

Size		Art.-No.
extra small	(5 - 5½)	01 395
small	(6 - 6½)	01 396
medium	(7 - 7½)	01 397
large	(8 - 8½)	01 398
extra large	(9 - 9½)	01 399

Packing: 150 pieces per dispenser box  
10 dispensers per outer carton



### BETATEX textured

Size		Art.-No.
extra small	(5 - 5½)	01 375
small	(6 - 6½)	01 376
medium	(7 - 7½)	01 377
large	(8 - 8½)	01 378
extra large	(9 - 9½)	01 379

Packing: 150 pieces per dispenser box  
10 dispensers per outer carton



## Chemotril

Extra strong, extra long Nitrile disposable gloves with extra mechanical strength and for high risk areas. Textured fingers and a beaded cuff add to extra resilience of this gloves.

Ambidextrous, comfortable fit, strong tear resistance, free of latex and latex proteins, powderfree. Resistant to almost all **cytotoxic chemicals**.

Free of silicone make the gloves suitable for many industrial applications. Excellent conductivity.



### CHEMOTRIL

Size		Art.-No.
small	6 - 6 ½	01 291
medium	7 - 7 ½	01 292
large	8 - 8 ½	01 293
x-large	9 - 9 ½	01 294
xx-large	10 - 10 ½	01 295

Packing: 50 pieces per dispenser box  
10 dispensers per outer carton





# Vinyl

Polyvinylchlorid (PVC), more commonly known as Vinyl, is a thermoplastic polymer. In order to make Vinyl usable for gloves and other products, it is mixed with plasticizers. It does not contain any Latex proteins. Vinyl is cheap to manufacture, but the decision for a Vinyl glove should be taken with care. Vinyl is not very elastic, which can negatively influence fit and resistance against biological hazards. Vinyl gloves should only be used in areas with a low risk factor. All Sanger disposable gloves made of Vinyl are ambidextrous, thiuram-free and have a beaded cuff edge.


## Pretex for light work



### PRETEX *smooth*

Size		Art.-No.
small	(6 - 6½)	01 126
medium	(7 - 7½)	01 127
large	(8 - 8½)	01 128
extra large	(9 - 9½)	01 129

Packing: 100 pieces per dispenser box  
10 dispensers per outer carton




The Pretex glove is made out of soft, tear-resistant Vinyl. It is suited for lighter work, like for example nursing care but can also be used in combination with chemicals.

- Latex Free
- Powder Free
- Smooth Surface



#### Standards:



# Polyethylene

Polyethylene (PE) is a versatile thermoplastic polymer. In its basic form it is transparent. PE gloves are manufactured by sealing two thin film PE foils to the shape of a glove. PE gloves are especially suitable for the food industry. Using PE gloves in areas with a high risk factor is not recommended. All disposable gloves made of polyethylene are ambidextrous and thiuram-free.

## Polyethylene Disposable Gloves


These economical gloves can be used for light work. The glove consists out of two heat-sealed, transparent polyethylene films with strong sealed seams. It does not influence the temperature perception.



### POLYETHYLENE

Size	Art.-No.
Women/29 cm	01 160
Men/29 cm	01 161

Packing: 100 pieces per polybag  
with dispenser opening  
100 polybags per outer carton



#### Standards:



## Pretex30 for protection from fluids

The Pretex30 glove is made out of soft, tear-resistant Vinyl. This makes it suitable for areas like nursing care, especially when in contact with liquids. It can be used in combination with chemicals.

- Latex Free
- Powder Free
- Smooth Surface


#### Standards:



### PRETEX30 *smooth*

Size		Art.-No.
small	(6 - 6½)	01 526
medium	(7 - 7½)	01 527
large	(8 - 8½)	01 528
extra large	(9 - 9½)	01 529

Packing: 100 pieces per dispenser box  
10 dispensers per outer carton



## Polyethylene Disposable Gloves *extra long, 55 cm*


This 55 cm long polyethylene glove reaches up to the elbow. That makes it very useful when protection for the entire forearm is needed.



### POLYETHYLENE *extra long*

Size/Length	Art.-No.
Universal 55 cm	01 165

Packing: 100 pieces per polybag  
with dispenser opening  
50 polybags per outer carton




#### Standards:



## Polyethylene Disposable Gloves *shoulder length, 94 cm*

### POLYETHYLENE *shoulder length*

Size/Length	Art.-No.
Universal 94 cm	01 164
Packing:	100 pieces per box 10 boxes per outer carton



This 94 cm long, soft and smooth polyethylene glove serves very well as protection for arm and clothes against soiling. The long glove is appropriate for veterinary use with extra long seams.

Standards:



## Veterinary Exam Gloves *with body protection, 95 cm*

The 95 cm long veterinary examination gloves have an integrated body protection to prevent the shoulder- and upper body section from soiling. One glove can be worn either on the left or the right side.

Sizes:  
Glove: Circumference at shoulder 76 cm x length 95 cm.  
Body: Width 60 cm (on either side) x length 78 cm.

With 2 strips for tying at front and back of approx. 70 cm length.



### POLYETHYLENE *with Body Protection, 95 cm*

Size/Length	Art.-No.
Universal 95 cm	03 165
Packing:	1 pieces per polybag 100 polybags per box 5 boxes per outer carton



Standards:



## Wall Shelf



An aluminium shelf for polyethylene gloves in women and men size. It can be adhered to smooth surfaces by double sided adhesive tape. Light-weight and complete with wall-mounting set.

### WALL SHELF

Art.-No.  
04 917

Size: Width 24 x height 15 x depth 1.5 cm

## Glove Dispenser Box Holders

- space saving for free standing at the place of work
- for assembly to wall or work top

Only minimal space is required for the box holders. The box holders allow a convenient use at work. All dispenser boxes of Sänger non-sterile disposable gloves fit into the white metal rack. Can be fixed on a surface or mounted on a wall with included set.



### DISPENSER BOX HOLDERS

Art.-No.	Capacity
04 911	1
04 912	2
04 913	3

### Dimensions:

1 tier L: 25.3 cm x H: 9.5 cm x D: 8.8 cm  
2 tiers L: 23.5 cm x H: 22.5 cm x D: 10 cm  
3 tiers L: 23.5 cm x H: 35.5 cm x D: 10 cm

## Hand Care Centre



The white metal holder for dispenser boxes of Sänger disposable gloves has additional space to store hand care products like OP-SeptBasic (page 27) or Sänger Derma Lotion (page 30). The Hand Care Centre allows easy access at work or at home. Can be positioned on a surface or mounted on a wall with included set.

### HAND CARE CENTRE complete with wall mounting set

Art.-No.  
04 914

Size of metal hand care centre:  
L: 37.5 cm x H: 8.5 cm x D: 9 cm

Size of metal hand care centre  
with dispenser box:  
L: 37.5 cm x H: 13 cm x D: 9 cm





## Duo Sterile for examinations

STERILE R

- Natural rubber latex
- Packed in pairs
- Polymer coating
- Powder free
- Textured surface
- Ambidextrous


DUO STERILE packed in pairs

Size	Art.-No.
small	01 206
medium	01 207
large	01 208
extra large	01 209

**Packing:**  
2 pieces of gloves packed in a medical grade paper wrapper and sealed in a sterile pouch.  
50 pairs per dispenser box  
4 dispensers per outer carton



Standards:




## Alfatex Sterile for examinations

- Nitrile
- Latex free
- Packed in pairs
- Powder free
- Textured surface
- Ambidextrous

ALFATEX STERILE packed in pairs

Size	Art.-No.
small (6 - 6½)	01 176
medium (7 - 7½)	01 177
large (8 - 8½)	01 178
extra large (9 - 9½)	01 179


**Packing:**  
2 pieces of gloves sealed in a sterile pouch.  
50 pairs per dispenser box  
10 dispensers per outer carton



Standards:



STERILE EO



## Pluto for surgery

The natural rubber latex with a polyurethane coating of the Pluto glove offers excellent tactile properties and guarantees a safe use with instruments, even in wet environment.

- Natural rubber latex
- Powder free
- Polyurethane coating
- Micro-textured surface
- Fully anatomically shaped
- Sterile

PLUTO powder free


Size	Art.-No.
6	01 310
6½	01 311
7	01 312
7½	01 313
8	01 314
8½	01 315

**Packing:**  
1 pair of gloves right/left packed in a medical grade paper wrapper and sealed in a sterile pouch.

50 pairs per dispenser box  
4 dispensers per outer carton



Standards:



STERILE R

## Eko for surgery

The Eko glove has the same material as the Pluto glove. The glove is lightly powdered with non-allergic food grade corn starch to provide easy donning.

- Natural rubber latex
- Lightly powdered
- Polyurethane coating
- Micro-textured surface
- Fully anatomically shaped
- Sterile

EKO lightly powdered

Size	Art.-No.
6	01 350
6½	01 351
7	01 352
7½	01 353
8	01 354
8½	01 355

**Packing:**  
1 pair of gloves packed in a medical grade paper wrapper and sealed in a sterile pouch.  
50 pairs per dispenser box  
4 dispensers per outer carton



Standards:





## Nuzone X2 for surgery (hypoallergenic)

The Nuzone X2 glove is ideal for sensitive skin. The synthetic polymer material offers excellent tactile properties, reduces hand perspiration and prevents allergies. These gloves can be used with cytostatics or chemicals. Free of MBT, Thiuram, Thiazol, and Thioureas.

- Synthetic Polymer
- Latex free
- Powder free
- Micro-textured surface
- Anatomically shaped
- Sterile



STERILE R

### NUZONE X2

Size	Art.-No.
6	01 250
6½	01 251
7	01 252
7½	01 253
8	01 254
8½	01 255

#### Packing:

1 pair of gloves right/left packed in a medical grade paper wrapper and sealed in a sterile pouch.

50 pairs per dispenser box  
6 dispensers per outer carton



#### Standards:



## Maxitex Duplex PF for high-risk surgeries

STERILE R

Duplex glove technology offers a primary yellow glove to be worn underneath a brown, non reflecting secondary glove. This ensures minimum pressure to the hands for long hours of continuous work without any signs of significant hand fatigue. They ensure highest safety without the outer glove slipping.



### Maxitex DUPLEX PF

Size to order	Secondary glove	Primary glove	Art.-No.
6	6	5½	01 610
6½	6½	6	01 611
7	7	6½	01 612
7½	7½	7	01 613
8	8	7½	01 614
8½	8½	8	01 615

#### Packing:

2 pairs of gloves packed in a medical grade paper wrapper and sealed in a sterile pouch.

25 x 2 pairs per dispenser box  
4 dispensers per outer carton



- Natural rubber latex 100%
- Electronically tested
- Traumatology and orthopaedic surgery
- High risk and reconstructive surgery
- Double shield protection
- Special PU surface coating

#### Standards:



## Maxitex Neuro PF for micro-surgery

The Maxitex Neuro PF sterile surgical glove is made of natural rubber latex and coated with Polyurethan, offering excellent tactile properties. The micro-textured surface allows safe use with instruments, even in wet environments. The brown colour reduces light reflections.




STERILE R

### Maxitex NEURO PF

Size	Art.-No.
6	01 510
6½	01 511
7	01 512
7½	01 513
8	01 514
8½	01 515

#### Packing:

1 pair of gloves right/left packed in a medical grade paper wrapper and sealed in a sterile pouch.



- Natural rubber latex
- Powder free
- Polyurethane coating
- Micro-textured surface
- Fully anatomically shaped
- Sterile

#### Maxitex NEURO Material Thickness:

Finger-tip: min. 0.17 mm  
Palm: min. 0.18 mm  
Cuff: min. 0.18 mm

#### Standards:



## Gynaecology for births


STERILE R

This extra long and extra thick sterile examination glove is suited for gynaecological areas, like underwater births. A special surface treatment makes it easy to don and gives very good tactile properties. The textured surface on hand and fingers provide good grip.



- Natural rubber latex
- Elbow length
- Powder free
- Structured surface
- Inner and outer chlorinated
- Anatomically shaped
- Sterile

#### Standards:



### GYNAECOLOGY

Size	Art.-No.
small (6 - 6½)	01 256
medium (7 - 7½)	01 257
large (8 - 8½)	01 258

#### Packing:

One pair of gloves packed in a medical grade paper wrapper and sealed in a sterile pouch.

25 pairs per dispenser box  
4 dispensers per outer carton





# Neutron 35

## NEUTRON 35

Size	Art.-No.
small	11 571
medium	11 572
large	11 573
x- large	11 574
xx-large	11 575

Packing:  
1 pair per polybag  
12 pairs per master polybag  
144 pairs per outer carton



Nitrile gloves for sterilizing Instruments. Highly resistant against many cleaning agents, chemicals and disinfectants.

- Nitrile
- Textured surface
- Anatomically shaped

Standards:



# MagnoGrip



Rubber gloves with special grip for easy donning of medical compression stockings. Protects against unwanted damages to the socks.

- Natural rubber latex
- Specially textured surface
- Anatomically shaped
- Inner flocklining

Standards:



## MAGNOGRIP

Size	Art.-No.
small	11 201
medium	11 202
large	11 203
extra large	11 204

Packing:  
1 pair per polybag  
10 pairs per master polybag  
50 pairs per outer carton



# Household Gloves

Rubber gloves suitable for working in the housekeeping area. A doubled coating grants extra protection and it can be used with a variety of disinfectants and chemicals.

- Natural Rubber Latex
- Textured Surface
- Anatomically Shaped
- Velour inner lining



Standards:



## PRIMAGRIP with zig-zag cuff

Size	Art.-No.
small	11 101
medium	11 102
large	11 103
extra large	11 104

Packing: Each pair is packed in a polybag with Euro-hole and EAN bar-code. 10 pairs per master polybag



# Radiation Protective Gloves

## Protech Proguard RR™

- Natural Rubber Latex
- Contains lead oxide
- Powder free
- Anatomically shaped
- Thickness 0.22 mm +/- 0.02 mm
- Textured Surface
- Re-sterilisable up to 4 times
- Very good tactile properties
- Tear resistant
- Straight cuff edge
- Low protein content, biocompatible

Extremely thin walled gloves with superior tactile properties as personal protection equipment against ionizing x-rays for use in pathology, orthopaedics, catheterizations, cardiology, intensive care, radiology, nuclear medicine, neuro-radiology, urology, endoscopy, resuscitation etc. According to directive 89/686/EEC.

Size, softness, elasticity, comfort and flexibility are on the level of surgical gloves. Four times re-sterilizable according to instruction for re-sterilization (1 instruction per dispenser box).


### Protech PROGUARD RR-1

Size	Art.-No.
6	01 710
6½	01 711
7	01 712
7½	01 713
8	01 714
8½	01 715


Packing:

One pair of gloves right/left packed in a medical grade paper wrapper and sealed in a sterile pouch.

5 pairs per dispenser box



Standards:



Reduction of Ionizing Radiation:

Type	60 KVP	80 KVP	100 KVP	120 KVP
RR-1	45%	35%	26%	23%
Radiaxon	59%	49%	42%	-

Request for full description and technical details.

Material Properties:

Type	Thickness	Lead Equivalence
RR-1	0.22 ± 0.02	0.016 mm Pb
Radiaxon min.	0.30	0.03 - 0.04 mm Pb

# Radiaxon

- natural rubber latex
- lead free and non-toxic
- powder free
- anatomically shaped
- thickness: min. 0.30 mm
- micro-textured surface
- for single use
- superior tactile properties
- tear resistant
- rolled cuff edge
- low protein contents

Sterile disposable glove mainly used for diagnostic heart catheterization, coronary angioplasty, angiocardiology, gynaecology, urology, orthopaedics etc.

The purpose of this radiation protective glove is to protect the hands from scattered secondary radiation exposure originating from the x-ray beam during fluoroscopic procedures. This glove is not intended for use in the direct or primary x-ray beam. Excellent grip on instruments. For single use.

## RADIAXON

Size	Art.-No.
6½	01 911
7	01 912
7½	01 913
8	01 914
8½	01 915
9	01 916

Packing:

1 pair of gloves right/left packed in a medical grade paper wrapper and sealed in a sterile pouch.

5 pairs per dispenser box



Standards:






## Säger ACTIVE Protect

For working safely with sharp instruments. Gloves have a safe and secure tight fit and reduce the danger of injury from sharp blades, edges or instruments drastically. Made of cut resistant yarn. Ideal for use in the sterilization department. Can be sterilized to be worn under medical gloves.

An inner label can be cut or marked with every washing cycle. Steam sterilization possible at 122.5°C after each washing cycle. Up to 4 times washable and sterilizable.



### AKTIVE PROTECT, non sterilized

Size	Art.-No.
medium (green)	01 585
large (grey)	01 586
extra large (yellow)	01 587

**Packing:** each pair in polybag  
30 pairs per carton

### AKTIVE PROTECT, sterilized by 122°

Size	Art.-No.
medium (green)	01 581
large (grey)	01 582
extra large (yellow)	01 583

**Packing:** each pair in sterile pouch  
30 pairs per carton

## Glove Liners



Seamless knitted nylon gloves, ambidextrous are comfortable to wear alone or inside other gloves. Reduced sweating underneath examination gloves. Tactile properties do not become affected due to half fingers and open, uncovered finger tips. Colour: white

- 100% Nylon
- Machine washable
- Ambidextrous

### NYLON KNITTED UNDERGLOVES

Size	Art.-No.
Universal	01 588

**Packing:** each pair in a polybag

## Cotton Knitting Gloves

Knitted cotton gloves with excellent fit. Particularly suitable for extra sensitive skin. Ideal for working with warm objects and protection against blunt edges. Good tactile properties as the sewing seams are situated at the side of the fingers.



### COTTON KNITTED GLOVES with fourchette fingers

Size	Art.-No.
small	01 415
medium	01 416
large	01 417
extra large	01 418
xx-large	01 419

**Packing:** each pair is packed in a polybag  
10 pairs per master polybag  
600 pairs per outer carton

### COTTON, ECONOMY ambidextrous

Size	Art.-No.
small	01 496
medium	01 497
large	01 498
extra large	01 499
xx-large	01 500

**Packing:** 12 pairs per polybag  
600 pairs per outer carton

### COTTON KNITTED GLOVES pairs with fourchette fingers

Size	Art.-No.
small	01 405
medium	01 406
large	01 407
extra large	01 408
xx-large	01 409

**Packing:** each pair is packed in a polybag  
10 pairs per master polybag  
600 pairs per outer carton

## Hand Scrub Brush

Reusable hand scrubbing brush, can be sterilized. Fits into standard wall holders.



### HAND SCRUB BRUSH sterilisable

Art.-No.
04 246

**Packing:** 30 pieces per dispenser box  
6 x 30 pieces per carton

## Disinfection Tray

Disinfection tray with lid and perforated inner plastic compartment. Additional lid to subdue light weight instruments.



### DISINFECTION TRAY

Capacity	Size	Art.-No.
0.1 Liter for small parts	dia Ø 85 x 40 mm	04 295
1 Liter	223 x 149 x 91 mm	04 291
3 Liter	300 x 210 x 115 mm	04 290
5 Liter	315 x 206 x 125 mm	04 292
10 Liter	496 x 328 x 195 mm	04 293
10 Liter for Endoscopes with outlet hose	870 x 235 x 160 mm	04 294



## OpSeptBasic

Ready-to-use alcoholic hand rub for hygienic and surgical hand disinfection, with excellent skin tolerance due to skin caring and moisturizing additives. Especially for medical departments and all working areas with an increased risk of infection and contamination. The skin tolerance is maintained even on long term and repeated usage. Dermatologically tested with result "excellent" 2013 Prof. Dr. med. Werner Voss.

**Spectrum of effect:** Bactericidal, fungicidal, tuberculocidal, virus-inactivating against Rota, Adeno, BVDV (HCV), Vaccinia and Norovirus. Also effective against MRSA: Effective within 15 seconds (tested acc. to EN 14476). Surgical hand disinfection: min 3 ml/90 seconds (acc. EN 12791). Observe user instructions.

Biozid register no. N-43561.



### OPSEPTBASIC

Size	Art.-No.
150 ml	04 277
500 ml	04 278
1 Liter	04 279
5 Liter	04 280





SangerDesinfekt


**Disinfectant for Instruments**

Is an aldehyde- and phenol-free, highly effective instrument disinfectant with excellent cleansing properties. An innovative combination of active substances guarantees a comprehensive, quick and material caring efficiency with very low usage concentrations. Is suitable for the reprocessing of all thermostable instruments and instruments, sensitive to temperature. Especially also for anaesthetic accessories, flexible endoscopes and endoscope accessories. Suitable for the manual reprocessing by immersion and for the ultra sonic bath as well as for the reprocessing of flexible endoscopes with the circulation method in so called half- and full automates.

Bactericidal, fungicidal (incl. *Aspergillus niger*), tuberculocidal (incl. *M. Terra*). Virus inactivating (covered viruses) BVDV/Vaccinia. HBV/HIV effective. Excellent material compatibility. Tested according to the standard methods of the DGHM (German Society of Microbiology and Hygiene), the EN norms for disinfection of instruments and according to the recommendations of the Robert Koch Institute for the efficiency against viruses. No odour.

**INSTRUMENT**

Size	Art.-No.
2 litres	04 275
10 litres	04 276



SangerDesinfekt


**SPRAY Disinfectant**

Instant disinfectant based on various alcohols and a new bispyridiamin derivate. For all areas with risk of infection and use in the vicinity of patients. Highly effective and extremely quick ready-to-use disinfectant in all clinical areas, laboratories, surgeries, dental areas, ambulances, for alcohol resistant surfaces and instruments. Perfect in-between surgical treatment and examinations when extremely short contact times are required. Bactericidal (incl. *Mycobacterium avium*)&TBC, fungicidal (incl. *Aspergillus niger*),

**SPRAY**

Size	Art.-No.
1 litres	04 282
10 litres	04 283

tuberculocidal (incl. *M-terrae*), virus-inactivating (incl. BVDV/HCV) Vaccinia incl HBV/HIV. According to DGHM standards. Effective against covered viruses according to current recommendations of RKI Robert Koch Institute, Germany. MRSA tested. Contact time: 60 seconds.



SangerDesinfekt


**SURFACE Disinfectant**

SangerDesinfekt SURFACE Disinfectant is an aldehyde- and phenol-free surface disinfectant based on a new quaternary ammonium combination with excellent cleansing characteristics. For all water- and wipe-resistant surfaces and for immersing objects. Ideal for areas with high degree of soiling and odour. Suitable also for kitchen and food-processing areas. Economical due to low dosage concentrations. Bactericidal, fungicidal, tuberculocidal (incl. *M.Terra*),

virus-inactivating (incl. HBV, HIV, BVDV-HCV, Vaccinia, Avian Flu virus). For disinfection of surfaces according to the standard methods of DGHM and current recommendations of RKI Robert Koch Institute, Germany.

**SURFACE**


Size	Art.-No.
2 litres	04 273
10 litres	04 274



SangerDesinfekt

**Disinfectant Wipes**

Hygienic ready-to-use wipes, soaked with Sanger-Desinfekt Spray Disinfectant for cleaning and disinfection of all alcohol resistant surfaces and equipments, also in the vicinity of patients. Economical due to the use of re-fill-boxes. The wipes are lint-free, none sticky, do not dry up and leave no visible residue. They are perfect for small



surfaces and small equipments. Effective against bacteria, fungi, tuberculosis (incl. *M.terrae*), inactivating viruses (incl. HBV/HIV). Do not contain Aldehyde. CE 0481.

**WIPES**

Size	Art.-No.
Refill set with 120 wipes	04 263
Box with wipes	04 265

**XXL Soft Wipes**

Extra large, well moisturizing wipes approx. 33 x 28 cm, dry. Just insert into the container and moisturize with approx. 1.2 – 1.5 litres of liquid disinfectant : i.e. SangerDesinfekt Spray, 10 litres, Art.-No. 04 283 (for information see page 28) or SprayOFF N as described below.

**Empty container** with reclosable opening. To be refilled with dry wipe rolls which can be moisturized with various liquid disinfectants.

**XXL SOFT WIPES**

	Art.-No.
70 Dry wipes Empty container	04 288
Empty container to insert fillings	04 289
Dry wipes: 1 roll per sachet, 8 sachets per carton	

**SprayOFF N**

Ready to use fast-acting disinfection solution for cleaning surfaces and medical equipment. Residue free and ideal for use between treatments. Effective against: BVDV (HCV): conc. / 30 sec., Vaccinia: conc. / 30 sec., MRSA effective and tested



(1 min), bactericide, fungicide, tuberculocide, virus-inactivating against all hulled virus, BVDV/(HCV) incl. HBV/HIV and Vaccinia according to recent suggestion by RKI, Noroviruses tested, 3 min.

**SPRAY-OFF N**

Size	Art.-No.
5 litres	04 285
10 litres	04 286

**Accessories****Description**

- 1 Spray bottle (empty)
- 2 Canister spanner
- 3 Canister pump
- 4 Canister tap

**ACCESSORIES**

No.	Art.-No.
1	04 903
2	04 902
3	04 904
4	04 901



### Derma Lotion

The ultimate skin protector against countless skin abusers and irritants: Contact allergies, harsh soaps, disinfectants, bleach, chlorine, salt water, acids, alcohol, unwanted odours, etc.

DERMA LOTION 250 ml
Art.-No.
02 231

### Aloe Vera Skin cleansing wipes

Aloe Vera skin cleansing wipes gently clean and moisturize the skin. Soft wet wipes are moisturized with a non irritating, alcohol free liquid containing Aloe Vera, Chamomile, Dead Sea minerals and Ginseng extract.



ALOE VERA WIPES	Art.-No.
Refill Pack (80 wipes)	04 244
Dimensions of wipe: approx. 20 x 16 cm	

### Lotus Washing Lotion

Lotus is an exclusive, soapless washing lotion. Keeping your skin soft and moisturized. Ideal for persons with sensitive skin.



LOTUS	
Size	Art.-No.
500 ml	04 240
1000 ml	04 241
10 litres	04 242

### Aloe 94% Gel

94% of this gel is pure Aloe Vera, making it one of the purest Aloe Vera skin care products in the world. Aloe Vera enhances skin regeneration. Contains approximately 200 bio-active substances including Vitamins B12, C & E, Beta-carotene, minerals, enzymes, sugars, anti-inflammatory salicylic acids and fatty acids. Recommended for pimples, minor burns, frostbite, insect bites, nappy rash, sunburns, dandruff etc. Also recommended as an effective soothing gel after shaving and hair removal.

For optimum results apply Aloe 94% several times daily on clean skin.



ALOE 94% GEL	
Contents	Art.-No.
100 ml	04 632

### Wall Mounted Lotion Dispenser



For washing lotion, liquid disinfectant and hand care lotions. Supplied with an empty dispenser. Suitable for Euro containers.  
**Material:** Body made of aluminium. Can be disinfected in an autoclav.  
**Lever:** Long stainless steel lever for operation by elbow.

LOTION DISPENSER		
Size	Art.-No.	Dimension
500 ml	04 920	82 x 283 x 82 mm
1000 ml	04 919	94 x 327 x 92 mm
Packing: 1 pieces		

### Noni

#### Soft and gentle skin care in just 20 minutes!

Cotton fibre layer inside the glove is moisturized with 60 ml liquid and coated with a thin polyethylene film, preventing moisture drying out.

Contents:  
**Noni extract (Morinda citrifolia):** cleans and re-vitalizes your skin.  
**Acacia Collagen:** with antiseptic and skin smoothing effect.  
**Ginkgo leaf extract:** provide relief of stressed skin, brings freshness and protects the skin keeping it healthy.



NONI SKIN CARE GLOVES
Art.-No.
04 410
NONI-FOOT CARE SOCKS
Art.-No.
04 411

**Macadamia nut oil:** is excellent for intensive skin care, rich in unsaturated fatty acids and vitamin A & E, high contents of palmitoleic acid which is a component of skin fat that promotes cell regeneration.

Packing: 1 pair per pouch






## Face Masks

### Printed

Sänger disposable surgical masks are made of a non-irritating, soft, light weight, non-woven material for optimum comfort. Ideal fit due to long, formable inner nose bridge. Sänger 3-layer surgical masks have a 99.6% filter effect and are free of glass particles and latex. Non-sterile.

Can also be worn comfortably while wearing glasses.  
Printed one-way caps see page 35.

Standards:  
 Med I

### FACE MASKS, PRINTED

#### with elastic ear loops

Size	Art.-No.
Child	02 137
Adult	02 136

Packing: 50 pieces per dispenser box  
6 dispensers per outer carton printed



## Face Masks

### Pink



### Standards:



#### WITH ELASTIC EAR LOOPS

Art.-No.  
02 103  
Packing: 50 pieces per dispenser box  
6 dispensers per outer carton pink

## Face Masks


### Blue

Sänger disposable surgical masks are made of a non-irritating, soft, light weight, non-woven material for optimum comfort. Ideal fit due to long, formable inner nose bridge. Sänger 3-layer surgical masks have a 99.6% filter effect and are free of glass particles and latex. Non-sterile.




#### WITH ELASTIC EAR LOOPS

Art.-No.  
02 108  
Packing: 50 pieces per dispenser box  
20 dispensers per outer carton blue

Standards:  
 Med I

#### WITH STRINGS FOR TYING

Art.-No.  
02 107  
Packing: 50 pieces per dispenser box  
20 dispensers per outer carton blue

Standards:  
 Med I



## Face Masks 3-layers

Sänger disposable surgical masks are made of a non-irritating, soft, light weight, non-woven material for optimum comfort. Ideal fit due to long, formable inner nose bridge. Sänger 3-layer surgical masks have a 99.6% filter effect and are free of glass particles and latex. Non-sterile. Can also be worn comfortably while wearing glasses.





#### WITH ELASTIC EAR LOOPS

Art.-No.  
02 102  
Packing: 100 pieces per dispenser box  
10 dispensers per outer carton green



#### WITH STRINGS FOR TYING

Art.-No.  
02 101  
Packing: 50 pieces per dispenser box  
20 dispensers per outer carton green


Standards:  
 Med I 



## Hygiene Masks



- with elastic loop
- with nose bridge
- non-glass
- suitable as secure dust filter mask

Standards:  
 Med I



### HYGIENE MASK


Colour	Art.-No.
white	02 106

Packing: 50 pieces per dispenser box  
20 dispensers per outer carton

## Paper Face Masks



- one layer
- with elastic ear loops
- 100 pieces per dispenser box with opening for easy removal

Standards:  
 Med I



### PAPER FACE MASKS


Colour	Art.-No.
white	02 105

Packing: 100 pieces per dispenser box  
100 dispensers per outer carton

## Mask Nik FFP3



Product features as mask Nik with additional comfort due to exhalation valve.  
EN 149:2001 A1:2009

Standards:  
 PSA I 

Nik Class	Art.-No.
FFP3	02 131

Packing: 10 pieces per box  
40 boxes per outer carton



## Sänger Face Masks

with eye shield and anti-glare-system

- Bacterial filter effect:  $\geq 98\%$
- Breathing resistance (Delta P)  $\leq 5.0$  mm
- Blood resistance 120 mm Hg acc. to ASTM F 1862
- Visor with anti-fog system, anti reflex edge, with 4 strings for tying

Can also be worn comfortably while wearing glasses.

Standards:



### WITH EYE SHIELD, ANTI-GLARE-SYSTEM

with strings for tying

Art.-No. 02 138

Packing: 25 pieces per dispenser box  
4 dispensers per outer carton

blue

Standards:



## Sänger Face Masks with eye shield

Sänger disposable surgical masks are made of a non-irritating, soft, light weight, non-woven material for optimum comfort. Ideal fit due to long, formable inner nose bridge. Sänger 3-layer surgical masks have a 99.6% filter effect and are free of glass particles and latex. Non-sterile.

Can also be worn comfortably while wearing glasses.



### WITH EYE SHIELD

with elastic ear loops

Art.-No. 02 104 132

Packing: 25 pieces per dispenser box  
4 dispensers per outer carton

blue

Standards:



### WITH EYE SHIELD

with elastic ear loops

Art.-No. 02 104 254

Packing: 25 pieces per dispenser box  
4 dispensers per outer carton

white

Standards:



## Surgical Caps

non-woven material

- air permeable
- excellent fit
- non-sterile
- size: medium: 50 cm / large: 54 cm circumference

Specially suitable as secure covering of hair in the medical trade. Hair and any other particles will be securely held back. The caps are soft and light weight. Air permeable material ensures reduction of perspiration even when wearing the caps for a long period. For one-way use.



### One-way Nurse Caps PLATE TYPE printed

One way nurse caps, white with blue print, with circular elastic.

Size	Art. -No.
universal	02 701

Packing: 100 pieces per polybag  
5 polybags per carton

### One-way Nurse Caps PLATE TYPE

One way nurse caps with circular elastic.

Size	Art.-No.
green	
medium	02 202
large	02 204
white	
medium	02 203
large	02 205
blue	
medium	02 206
large	02 207
pink	
medium	02 709
large	02 707

Packing: 100 pieces per polybag  
5 polybags per carton

### One-Way Nurse Caps PLEATED

Pleated nurse caps packed in polybags, easy removal. Air permeable material.

Size	Art.-No.
green	
medium	02 201
white	
medium	02 208
blue	
medium	02 209
pink	
medium	02 703
red	
medium	02 702

Packing: 100 pieces per polybag  
5 polybags per carton





## PE-Apron with Arm Sleeves

Strong water resistant polyethylene foil. Blue, textured surface, non-transparent. With neckband and two strings for tying at the waist. With small thumb hole to ensure gloves will cover the arm sleeves to securely keep out any splashing liquids or other contaminating material. Ideally suited for disinfecting instruments in liquid. For one-way use.

Standards:



### PE-APRON with arm sleeve

Dimension	Art.-No.
120 x 170 cm	02 251

**Packing:**  
Each piece in a polybag  
10 polybags in a master polybag  
10 master polybags per outer carton

## Disposable Gown

White, non-sterile, non-woven material for single use, i.e. visitors, medical staff or patients, etc. With two ties at the hip and at the neck to close the gown at the back. Each piece in a polybag.



### DISPOSABLE GOWN

Size	Art.-No.
x-large white	02 248
xx-large white	02 249

**Packing:**  
Each piece in a polybag  
50 pieces per outer carton



## Disposable Apron

Ideal zic-zac layers: Remove one apron and the next one will follow through the dispenser opening for next removal.

Size	Art.-No.
70 x 120 cm	02 307

**Packing:** 100 pieces per dispenser box  
10 boxes per outer carton

## Sleeve Covers

Made of strong polyethylene foil, with elastic at both ends.



### SLEEVE COVERS

Colour	Art.-No.
white	02 218
blue	02 219

**Packing:**  
50 pieces per polybag  
40 polybags per outer carton

## Shoe Cover Box

For donning shoe covers quickly. Easy to use and no contact with the contaminated area. Mechanical functionality, no power-source required, making it quick to set up. Ideal for areas with high hygienic demands, like hospitals, clean rooms, laboratories, etc. Refill pack contains 50 covers.

**Size:** approx. 47 x 25 x 19 cm  
**Weight:** 5.85 kg  
**Capacity:** 50 pcs.  
**Machine:** Silver colour metal plated plastic case with carrying handle.



SHOE COVER BOX, EMPTY  
Art.-No. 04 890



REFILL PACK SHOE COVERS  
Art.-No. 04 891

## Disposable Shoe Covers

1) Material: Non-woven. Strong and durable non-woven material, textured, with all seams sewn. Colour white. Non-sterile.

2) Material: PE foil. Textured surface tear resistant, water proof, colour blue. Non-sterile.  
Size: medium.

2) Material: CPE foil. Water resistant foil with textured slip resistant surface. Very durable, colour blue. Non-sterile.  
Size: large.

3) Material: Non-woven. Colour white. Disposable boot covers with textured surfaces. 40 cm height. Non-sterile.



### ONE-WAY SHOE COVERS

Material	Colour	Art.-No.
1) Non-woven	white	02 220
2) PE	blue	02 230
2) CPE	blue	02 234
3) Non-woven	white	02 232

**Packing:**  
100 pieces per polybag  
10 polybags per outer carton

## Colonoscopy Pants

Colonoscopy pants, non-sterile, with rubber band in hip area. With rear opening. Flap is held open or closed by Velcro® seal fasteners.

**Material:** Non-woven, non-transparent



### COLONOSCOPY PANTS

Size	Art.-No.
Blue, with rear slot	
large	02 123

**Packing:**  
10 pieces per polybag  
10 polybags per outer carton

## Disposable Panties

Made of non-woven material with elastic on waist and legs. Light, soft and air permeable. For one-way use. Reinforced section up to the waist. Non-sterile.



### DISPOSABLE PANTIES

Size	Art.-No.
XL (59 x 70 cm)	02 215
XXXL (66 x 75 cm)	02 216

**Packing:**  
50 pieces per polybag  
20 polybags per outer carton

## Elastic Mesh Panties

Made of elastic mesh material. Suited for mobile or bed-ridden patients. Holds firm and steady. One-size fits all with seams on the side to enlarge the size by tearing. Keeps stock quantities down.



### ELASTIC MESH PANTIES

Size	Art.-No.
universal	02 125

**Packing:**  
50 pieces per dispenser box  
500 pieces per outer carton

## Beard Cover

Non-woven, white, with rubber band, breathable.



### BEARD COVER

Size	Art.-No.
53 x 22 cm	02 237

**Packing:** 100 pieces in a polybag  
1000 pieces per outer carton

## Disposable Gown

White, non-sterile, non-woven material for single use, i.e. visitors, medical staff or patients, etc. With two ties at the hip and at the neck to close the gown at the back. Each piece in a polybag.

### DISPOSABLE GOWN

Size	Art.-No.
x-large	02 248
xx-large	02 249

**Packing:**  
Each piece in a polybag  
50 pieces per outer carton



## Disposable Overall Cosmo

Water resistant overall made of waterproof-coated polypropylene. Comfortable to wear. Several sizes for excellent fit. Available with elastics at cuff, hood and legs. Zip fastener at the front. Available in 4 sizes for comfortable fit. Soft, strong, non-fluff white colour material with very good resistance to abrasion. Protects against many chemicals, paints, dust etc. Ideal to have in any emergency body protection kit. Suitable for laboratories, industrial work.

CE

### DISPOSABLE OVERALL

Type	Art.-No.
small	02 110
medium	02 111
large	02 112
xx-large	02 115

**Packing:** Each piece in a polybag  
25 pieces per outer carton



## PE Gown

White colour transparent polyethylene gown with push buttons in front. Available with hood or with collar. Ideal protective clothing for visitors, for protection of cloth against wet, dirt and odour. Colour: white/transparent.

CE

### PE GOWN

Type	Art.-No.
PE gown with collar	02 113
PE gown with hood	02 114

**Packing:**  
PE gown: Each piece in a polybag  
10 pieces per master polybag  
500 pieces per outer carton

## Disposable Working Gown

Made of polypropylene, 35g/qm. Suitable for labs, visitors etc. With push buttons in front, collar and breast pocket. Colour white.



### WORKING GOWN

Size	Art.-No.
x-large	02 192
xx-large	02 193

**Packing:**  
5 pieces in a polybag  
50 polybags per outer carton







## Doctors Jacket, non-sterile

Disposable jacket for one-way use. Knitted cuffs and collar section for added comfort. Suitable as disposable jacket to cover clothing in the doctors daily practice, as warming jacket for anesthetists etc. Keeps private clothes clean. With push buttons in front.

Material: SMS Length: approx. 70 cm

### DOCTORS JACKET

Size	Art.-No.
Universal	02 247

Packing: Each piece in a polybag, 50 pieces per carton

## Easy & Clean Sticky Mats

Sticky Mats provide a unique, cost-effective solution to trap dirt before it is tracked into your workplace, office or clinic. Easy-to-use and long lasting. 30 adhesive sheets are individually removable from the pad.


Used mats can be disposed easily and quickly. The mats are numbered so replacements can be obtained in advance.

The underside is adhesive to keep the mat securely in place whilst in use.



### EASY & CLEAN

	Colour: Blue	Colour: Transparent
<b>Dimensions</b>	<b>Art.-No.</b>	<b>Art.-No.</b>
45 x 90 cm	04 887	04 897
60 x 90 cm	04 888	04 898
66 x 114 cm	04 889	04 899



## Doctors Gown, non-sterile

Disposable gowns for one-way use, with knitted cuff, velcro seal at the back. Additionally 4 strings for tying the gown at the side and at the back. Can be supplied in three below mentioned qualities. With and without wet protection zone in front and inside arms. All seams are sewn.

### Doctors Gown ECO non-sterile

Material: SB35 hydrophobic polypropylene non-woven textile, spunbond type, green.

without wet protection zone		with wet protection zone	
Size	Art.-No.	Size	Art.-No.
medium	02 480	medium	02 485
large	02 481	large	02 486
extra large	02 482	extra large	02 487
xx-large	02 483	xx-large	02 488

#### Packing:

Without wet protection zone: Each piece in a polybag, 50 pieces per carton  
With wet protection zone: Each piece in a polybag, 40 pieces per carton

### Doctors Gown STANDARD non-sterile

Material: SMS hydrophobic material made of 100% polypropylene non-woven textile 50 g/sqm, 3-layer structure spunbond-meltblown-spunbond, blue.

without wet protection zone		with wet protection zone	
Size	Art.-No.	Size	Art.-No.
medium	02 221	medium	02 227
large	02 225	large	02 228
extra large	02 226	extra large	02 229
xx-large	02 256	xx-large	02 259

Packing: Each piece in a polybag, 40 pieces per carton

#### Length and sizes

Size	Length	Width
medium	125 cm	150 cm
large	135 cm	155 cm
extra large	145 cm	164 cm
xx-large	155 cm	164 cm



## Sänger Surgical Gown, sterile

### One-way Surgical Gown ECO sterile

For short, simple procedures or surgical assistants, green. Material: SB35

without wet protection zone		with wet protection zone	
Size	Art.-No.	Size	Art.-No.
medium	02 490	medium	02 495
large	02 491	large	02 496
extra large	02 492	extra large	02 497
xx-large	02 493	xx-large	02 498


Packing: Each piece in a sterile pouch with wrapping sheet  
without wet protection zone: 30 pieces per carton  
with wet protection zone: 30 pieces per carton

### One-way Surgical Gown STANDARD sterile

For surgical procedures of all kinds, blue. Material: SMS

without wet protection zone		with wet protection zone	
Size	Art.-No.	Size	Art.-No.
medium	02 240	medium	02 243
large	02 241	large	02 244
extra large	02 242	extra large	02 245
xx-large	02 253	xx-large	02 255

Packing: Each piece in a sterile pouch with wrapping sheet  
without wet protection zone: 28 pieces per carton  
with wet protection zone: 25 pieces per carton




#### Available materials:

**SB35**  
Hydrophobic polypropylene non-woven textile, spunbond type. For low fluid procedures.

**SMS**  
Hydrophobic material made of 100% polypropylene non-woven textile 50 g/sgm, three-layer structure spunbond-meltblown-spunbond. For low fluid procedures.sweating.

## Bill Cover Sterile Surgical Retractor Covers



**Cotton surgical retractor hook covers for the safe, slip-free use of surgical retractors of all sizes.**

- X-ray identification marking
- Secure attachment of the cover to retractor hook
- For single use
- Coiled, in individual sterile packing
- Ready for use and latex-free
- Individual sets available
- Individual sizes available

Specially good grip due to the double layer of the material on the retractor hook to prevent slipping on the bare metal. Suitable for liver hooks, lung spatulas or bowel clamps, etc. The covers are supplied sterile, ready for use. They are in accordance with the European Directive relating to medical products 93/42/EEC.


Free testsamples available upon request.



### Packing:

- |                         |                        |
|-------------------------|------------------------|
| <b>1 pair:</b>          | <b>1 set:</b>          |
| 2 retractor hook covers | 1 retractor hook cover |
| 2 adhesive strips       | 2 adhesive strips      |
| 50 pairs per carton     | 50 sets per carton     |

Size	Dimension	Art.-No.	Unit
0	10 x 55 mm	02 640	1 Pair
1	10 x 120 mm	02 644	1 Set
21	25 x 150 mm	02 645	1 Set
22	25 x 200 mm	02 646	1 Set
2	25 x 250 mm	02 641	1 Set
23	25 x 300 mm	02 648	1 Set
7	25 x 350 mm	02 647	1 Set
3	25 x 400 mm	02 642	1 Set
4	25 x 500 mm	02 650	1 Set
6	35 x 350 mm	02 651	1 Set



STERILE EO

Art.-No.	02 350
Size	50 x 50 cm
Aperture	oval 6 x 8 cm
Pcs./Pack	1

Art.-No.	02 352
Size	50 x 70 cm
Aperture	round 10 cm diameter
Pcs./Pack	1

Art.-No.	02 354
Size	75 x 90 cm
Aperture	round 10 cm diameter
Pcs./Pack	1

Art.-No.	02 356
Size	75 x 120 cm
Aperture	round 10 cm diameter
Pcs./Pack	1

Art.-No.	02 357
Size	100 x 150 cm
Aperture	oval 6 x 8 cm
Pcs./Pack	1

## Surgical Drapes, sterile For Patients and Inventory

Material: hydrophobic SMS polypropylene nonwoven material 50 g/sqm, 3-layered material, type: meltblow-spunbond-meltblown

- for low fluid procedures

### SURGICAL DRAPES Material: SMS

Art.-No.	Size (Drape)	Adhesive tape	Pcs./Pack
02 270	50 x 50 cm	-	1
02 271	50 x 70 cm	-	1
02 273	75 x 75 cm	-	1
02 274	75 x 90 cm	-	1
02 283	75 x 75 cm	2.5 x 75 cm	1
02 286	100 x 120 cm	5 x 80 cm	1
02 287	150 x 180 cm	8 x 80 cm	1


## Sänger Surgical Barrier System

It is a system of mutually compatible elements, enabling effective prevention against infection during surgical procedures. It safely isolates any possible source of infection, especially the patient's body, the operating theatre staff and equipment, and simultaneously absorbs body fluids, which may carry any infection.

The system contains

- surgical drapes
- adhesive tapes
- complementary drapes
- table covers, MAYO covers
- incision films
- fluid collection pouches
- sleeve covers
- covers for device

Request for our complete catalogue on surgical drapes and surgical disposables.



### Available Materials


- **DS2**  
Two layer laminated material  
Upper layer: absorbent viscose 28 g/sqm  
Under layer: impervious polyethylene foil 25 µm, for medium fluid procedure.
- **DS3**  
Three layer laminated material  
Upper layer: absorbent viscose 28 g/sqm  
Intermediate layer: impervious polyethylene foil 25 µm. Under layer: absorbent non-woven textile, for high fluid procedure
- **SMS**  
Hydrophobic material made of 100% polypropylene non-woven textile 50 g/sqm, three layer structure spunbond-meltblown-spunbond, for low fluid procedure

### SURGICAL DRAPES Material: PF2

Art.-No.	Size (Drape)	Adhesive tape	Pcs./Pack
02 290	50 x 50 cm	-	1
02 291	50 x 70 cm	-	1
02 294	75 x 90 cm	-	1
02 301	50 x 50 cm	2.5 x 50 cm	1
02 302	50 x 70 cm	2.5 x 50 cm	1
02 305	75 x 90 cm	2.5 x 90 cm	1

- **SBPE**  
Two layer laminate  
Upper layer: polypropylene non-woven textile 30 g/sqm  
Under layer: impervious polyethylene foil 25 µm
- **PF2**  
Two layer laminate  
Upper layer: absorbent tissue 20 g/sqm  
Under layer: impervious polyethylene foil 25 µm
- **SB**  
Hydrophobic polypropylene non-woven textile, spunbond type
- **PE**  
Polyethylene foil

Depending on the type and kind of operation all articles can be supplied in different types of materials.




## Surgical Drapes, sterile For Patients and Inventory

Material: PF2, 2-layered laminated material, liquid resistant  
Upper layer: absorbent cellulose 20 g/sqm  
Under layer: impervious polyethylene foil 20 µm

- for low fluid procedures



Fenestral Drapes with adhesive area Material DS2




Art.-No.	02 361
Size	45 x 75 cm
Aperture	oval 6 x 8 cm
Pcs./Pack	1

Art.-No.	02 365
Size	75 x 90 cm
Aperture	round 10 cm diameter
Pcs./Pack	1

Art.-No.	02 367
Size	90 x 110 cm
Aperture	round 7 cm diameter
Pcs./Pack	1

Fenestral Drapes Material DS3




Art.-No.	02 370
Size	50 x 70 cm
Aperture	oval 6 x 8 cm
Pcs./Pack	1

Art.-No.	02 372
Size	90 x 120 cm
Aperture	oval 6 x 8 cm
Pcs./Pack	1

Art.-No.	02 373
Size	120 x 150 cm
Aperture	oval round 8 cm diameter
Pcs./Pack	1

Surgical Drapes Material DS2



Art.-No.	02 310
Size drape	50 x 50 cm
Pcs./Pack	1

Art.-No.	02 311
Size drape	50 x 70 cm
Pcs./Pack	1

Art.-No.	02 312
Size drape	40 x 60 cm
Pcs./Pack	1

Tube cover polyethylene




Art.-No.	02 475
Size	8 cm x 120 cm
VE	1

Adhesive Tape

Self-adhesive surgical tape with absorbent surface area.



Material:	NV28
Art.-No.	02 457
Size	5 cm x 30 cm
Pcs./Pack	1



Incision Foil

Self adhesive incision foil for affixing the surgical drape onto the patients skin. On carrier paper for easy affixation at the correct place. Sterile packing. Material: soft, tear-resistant polyurethane, transparent.

Size	Art.-Nr.
40 x 50 cm	02 609

Custom packs:

We specialise in assembly of individual packs according to customers requirements, offering selected standard items for everyday use. Any item of this catalogue can be used for configuration of the custom pack considering: surgical techniques, organisation and equipment of operating theatres and number of members of the surgical team.

Surgical Drapes Material DS3


Art.-No.	02 331
Size drape	50 x 70 cm
Pcs./Pack	1

Art.-No.	02 334
Size drape	75 x 90 cm
Pcs./Pack	1

Art.-No.	02 336
Size drape	100 x 120 cm
Pcs./Pack	1

Ultrasonic Probe Protective Covers

Intended use: For vaginal and rectal probes  
Surface: smooth, dry, without lubricant.  
Can be delivered in various sizes, loose in a dispenser box or individually sealed.



LATEX, 144 pieces		
ø	Length	Art.-No.
30 mm	195 mm	02 666
Packing: 144 pieces per box		

LATEX		
ø	Length	Art.-No.
13.5 mm	55 mm	02 660
20 mm	200 mm	02 661
28 mm	200 mm	02 664
30 mm	300 mm	02 665
Packing: 100 pieces per polybag		


Ultra Sonic Gel



A proven contact gel for the exterior use of ultrasonics in therapeutical and diagnostical medicine. The gel is skin and mucus compatible, pH neutral and not sensitizing. Drip and bubble free.

ULTRA SONIC GEL

Size	Art.-No.
250 ml bottle	04 210
500 ml bottle	04 211
1 Liter bottle	04 212
5 Liter container	04 213
5 Liter cubitainer	04 214



Baby Blanket

- Natural Fiber
- Bio-degradable
- Pleasant to the skin
- Soft and smooth

No coating, therefore not slippery in wet or dry condition. Sterilized. Absorbs fluids and breathable.



BABY-BLANKET disposable

Size	Art.-No.
100 x 100 cm	04 506 155
Packing: Single units sterilized	

## Washable pad

Liquid-proof + absorbent

For adults and children: Four-layer reusable system

Machine washable underpad with soft surface. Will not crinkle when used to avoid pressure points on human skin, which cause wounds. Long lasting quilting seams. Can be boil washed repeatedly with short drying period in household or professional dryer. Washing instructions are provided. For use with elderly people, hospital patients, children. Will protect mattresses, chairs and furniture from soiling.


**Material:** Polyester / Viscose



Size:	Art.-No.
90 cm x 75 cm	03 780
90 cm x 85 cm	03 781

**Packing:**  
1 piece per polybag  
20 pieces per carton

**Absorption capacity:** approx. 1000 ml - 1200 ml



Size: 90 x 85 cm,  
cotton side flaps: approx. 42 cm



Size: 90 cm x 75 cm

## Mattress protection

Hygienic cover of mattresses and treatment tables. Rapid change and treatment possible, no washing required.  
Material: Polypropylene

### MATTRESS PROTECTION

Size:	Art.-No.
90 x 200 x 20 cm	04 511

**Packing:**  
60 pieces per carton



STERILE EO

## Disposable Needles

- suitable for cylindric ampule and Luer syringes
- with plastic hub, precision safety glued
- siliconized. For single use
- three edge lancet-grinding for painless injection
- particularly voluminous flow
- 30 G suitable for intraligamental anesthesia
- easy mounting of the needles due to precisely cut thread
- clean puncture into the septum without punching of rubber particles due to precision grinding



### DISPOSABLE NEEDLES


Diameter	Length	Colour	Art.-No.
<b>Injection</b>			
G 25/0.5 mm	42 mm	red	03 700
G 27/0.4 mm	42 mm	yellow	03 703
	37 mm	yellow	03 704
	23 mm	yellow	03 705
G 30/0.3 mm	23 mm	green	03 706
	13 mm	green	03 707
	10 mm	green	03 708

### Injection Luer

G 27/0.4 mm	42 mm	grey	03 712
	23 mm	grey	03 713

### Packing:

Each piece sterile packed  
100 pieces per box



### BLUNT NEEDLES


Diameter	Length	Colour	Art.-No.
<b>Metric Blunt</b>			
G 20/0.9 mm	32 mm	red	03 720
G 21/0.8 mm	23 mm	green	03 721

### Luer Blunt

G 20/0.9 mm	42 mm	red	03 722
	23 mm	yellow	03 723

### Packing:

Each piece sterile packed  
100 pieces per box



## Disposable Needles Blunt

- with plastic hub, precision safety glued
- suitable for cylinder ampule injections
- rounded tip with rounded edges for use in perio-dontal pockets (blunt needles) and in the root canal (blunt needles Luer Endo)
- easy mounting of the needles due to precisely cut thread
- clean puncture into the septum without punching of rubber particles due to precision grinding.
- Transcodent: CE 0482

### ENDO IRRIGATION NEEDLES

Diameter	Length	Colour	Art.-No.
<b>Irrigation Luer Endo</b>			
G 27/0.4 mm	25 mm	grey	03 732
G 30/0.3 mm	25 mm	yellow	03 733

### Packing:


Each piece sterile packed  
25 pieces per box

### Luer Endo

G 25/0.5 mm	42 mm	yellow	03 724
G 27/0.4 mm	23 mm	blue	03 725
G 30/0.3 mm	23 mm	green	03 726

### Packing:

Each piece sterile packed  
100 pieces per box



## Endo Irrigation Needles

Luer

- only irrigation needle with 2 side openings
- for effective irrigation of the root canal
- openings and rounded end prevent irrigation beyond the apex
- irrigation of the entire root canal without danger
- Luer or Luer lock single, blister packed
- for single use
- Transcodent: CE 0482






STERILE R

## Disposable Syringes

Luer (Slip) two part disposable syringe without needle. 2 ml syringe with centred tip, 5-20 ml offset tip. Transparent barrel with easy-to-read accurate calibration lines. Non-toxic material, free from latex, PVC and pyrogen. Single blister packed. Sterile.




### LUER (SLIP)

Volume	Art.-No.
2 ml	03 102
5 ml	03 105
10 ml	03 110
20 ml	03 120

**Packing:** 100 pieces/box  
 2 ml: 2,000 pieces/outer carton  
 5 ml / 10 ml: 1,200 pieces/outer carton  
 20 ml: 700 pieces/outer carton

## Disposable Syringes

Luer (Lock) tip. Two part disposable syringe. Without needle. Manufactured to ensure the perfect mating of barrel and plunger. Non-toxic material, free from latex, PVC and pyrogen. Single blister packed. Sterile.





### LUER (LOCK)

Volume	Art.-No.
2 ml	03 202
5 ml	03 205
10 ml	03 210
20 ml	03 220

**Packing:** 2 ml: 100 pcs./box, 2,800 pcs./outer carton  
 5 ml: 100 pcs./box, 1,900 pcs./outer carton  
 10 ml: 100 pcs./box, 1,200 pcs./outer carton  
 20 ml: 50 pcs./box, 750 pcs./outer carton

## Disposable Needles

Sterile disposable needles, designed for Luer (Slip) connection. Latex free, PVC free and sterile. Surface treatment and coating for easy permeation. The hub of the needle is made of transparent material, which facilitates visual control of the flow of medicine or blood through the needle. Thin walled for smooth flow of fluid. Each piece is in a sterile blister packed.



### DISPOSABLE NEEDLES

Gauge	Size	Art.-No.
26G x 5/8"	0.45 x 16 mm	03 416 (brown)
25G x 4/5"	0.50 x 20 mm	03 520 (orange)
25G x 1"	0.50 x 25 mm	03 525 (orange)
23G x 1"	0.60 x 25 mm	03 625 (blue)
22G x 1 1/5"	0.70 x 30 mm	03 735 (black)
21G x 1 1/2"	0.80 x 40 mm	03 840 (green)
20G x 1 1/2"	0.90 x 40 mm	03 940 (yellow)
19G x 1 1/2"	1.10 x 40 mm	03 140 (ivory)
18G x 1 1/2"	1.20 x 40 mm	03 1240 (pink)

**Packing:** 100 pieces per carton  
 4.000 pieces per outer carton



**Volume** 0.25 Litre  
**Height:** 130 mm  
**Lid:** 80 mm x 45 mm  
**Packing:** 200 pcs./carton  
 4.000 pcs./pallet



**Volume** 0.5 Litre  
**Height:** 157 mm  
**Lid dia.:** 86 mm  
**Packing:** 100 pcs./carton  
 2.000 pcs./pallet



**Volumen** 3,5 Liter  
**Height:** 314 mm  
**Lid dia Ø** 145 mm  
**Packing:** 50 pcs./carton  
 600 pcs./pallet



**Volume** 6 Litre  
**Height:** 240 mm  
**Width:** 280 mm  
**Depth:** 190 mm  
**Packing:** 24 pcs./carton  
 288 pcs./pallet

Available in different kinds of sizes and shapes. Suggested for syringes, scalpels, etc. for safe disposal. With different kinds of lids. Securely lockable for transportation and disposal. Comply to UN and ADR standards.

Autoclavable, sticker-free.

**For fixtures to table tops, brackets or trolleys: please ask for separate leaflet.**



**Volume** 1.5 Litre  
**Height:** 152 mm  
**Lid dia.:** 145 mm  
**Packing:** 90 pcs./carton  
 1.080 pcs./pallet



**Volumen** 4,0 Liter  
**Height:** 150 mm  
**Width:** 280 mm  
**Depth:** 170 mm  
**Packing:** 40 pcs./carton  
 480 pcs./pallet



**Volume** 11 Litre  
**Height:** 330 mm  
**Width:** 280 mm  
**Depth:** 190 mm  
**Packing:** 28 pcs./carton  
 336 pcs./pallet

## WASTE REMOVAL SAFETY CONTAINERS

Volume	Art.-No.	
0.25 Litre	04 926	pocket-size
0.5 Litre	04 921	standard lid
1.5 Litre	04 922	standard lid
2.2 Litre	04 918	standard lid
3.5 Litre	04 923	standard lid
4.0 Litre	04 929	standard lid
6.0 Litre	04 927	standard lid
6.0 Litre	04 924	with flap
11 Litre	04 928	standard lid
11 Litre	04 925	with flap



**Volume** 2.2 Litre  
**Height:** 204 mm  
**Lid dia.:** 145 mm  
**Packing:** 70 pcs./carton  
 840 pcs./pallet



**Volume** 6 Liter  
**Height:** 210 mm  
**Width:** 280 mm  
**Depth:** 190 mm  
**Packing:** 38 pcs./carton  
 456 pcs./pallet



**Volume** 11 Litre  
**Height:** 360 mm  
**Width:** 280 mm  
**Depth:** 190 mm  
**Packing:** 20 pcs./carton  
 240 pcs./pallet

## SängerPlast Cotton Roll Plaster

Cotton tape, strong fixation, air permeable. Skin coloured fixing tape. 100% cotton. SängerPlast Cotton Roll Plaster is tearable in length and width, without a scissor. Radiotransparent. Non-sterile.



### COTTON PLASTER

Size	Art.-No.	Rolls per box	Boxes per outer carton
1.25 cm x 9.14 m	03 532	24	10
2.5 cm x 9.14 m	03 533	12	10
5 cm x 9.14 m	03 534	6	10

## SängerPlast Pore

Air and steam permeable tape, hypoallergenic for delicate skin i.e. children and elderly, neck, face etc. Easy removal. Micro porous structure allows the skin to breathe. Can be easily torn by hand. Radiotransparent. Non-sterile.



### PORE

Size	Art.-No.	Rolls per box	Boxes per outer carton
1.25 cm x 9.14 m	03 522	24	10
2.5 cm x 9.14 m	03 523	12	10
5 cm x 9.14 m	03 524	6	10

## SängerPlast Minis

High performance first aid strips. Skin colour, skin friendly, air permeable. Cotton, for quick and reliable application. 250 pieces loose packed per box with lid. Specially suitable for blood sampling, injections, small wounds etc. Non-sterile.



### MINIS

Size	Art.-No.
20 x 40 mm	03 489
<b>Packing:</b>	
250 pieces per box	
20 boxes per outer carton	

## SängerPlast Maxis

High performance first aid strips. Skin colour, individually wrapped. Hypoallergenic adhesive avoids any skin irritation. 200 individually packed waterproof plaster strips, non-sterile in a box with lid. Pad will not adhere to the wound surface. Non-sterile.



### MAXIS

Art.-No.	
03 402	
<b>Contents:</b>	
20 x 70 mm	150 pcs.
20 x 60 mm	30 pcs.
20 x 40 mm	20 pcs.
<b>Packing:</b>	
200 pieces per box	
51 boxes per outer carton	



## SängerPlast Sensitive

Sensitive, for sensitive skin, painless removing

Hypoallergenic wound plaster for delicate skin sections such as face skin, neck, skin of children and elderly etc. Good adhesive, air permeable and no skin irritations. Air permeable through cell structure of the material. Easy removing. Pad does not stick to the wound surface. Hypoallergenic adhesive avoids any skin irritation. Radiotransparent. Roll can be cut to individual sizes. Non-sterile.



### SENSITIVE

Size	Art.-No.
6 cm x 5 m	03 465
<b>Packing:</b>	
1 roll per box	
24 boxes per outer carton	

## Sport-Tape

For non elastic tape bandages to stabilise joints and prevent sports injuries from all kinds of exercises. The tape can be torn without scissors lengthwise and across. Adheres well to skin and other surfaces and does not easily come off, even under sweaty or wet conditions. Does not leave trace particles behind when being removed from skin or from other surfaces i.e. guards, socks, shoes, etc. CE Made of 100% cotton with zinc-oxide-adhesive. Tested for anti-allergenic and non-irritational properties. Water resistant. Ideal for use as functional supporting bandage. Prevents sprains and protects from injuries during recovery period. Can be used in the treatment of luxations and is used as protective measure from injuries during sports activities. For Pre-Tape refer to page 53.



### SPORT-TAPE

Colour	Size	Art.-No.	Rolls/box	Rolls/carton
white	2.5 cm x 10 m	03 481	18	432
white	3.8 cm x 10 m	03 480	12	288
white	5 cm x 10 m	03 560	6	60
blue	2.5 cm x 10 m	03 483	18	432
blue	3.75 cm x 10 m	03 482	122	88
yellow	3.75 cm x 10 m	03 501	12	288
skin	3.75 cm x 10 m	03 503	12	288
red	3.75 cm x 10 m	03 505	12	288

## Lister Scissors for Bandages

First aid scissors according to DIN 58 279 suitable for hospitals, doctors practises, first aid services, emergency sets and boxes etc. The classical bandage scissor is made of stainless steel for professional use.



### LISTER SCISSORS

Art.-No.	
03 621	
<b>Packing:</b>	
200 pieces per outer carton	



# SängerPlast Classic Cotton

- highly absorbent
- anti-allergenic
- effective dirt barrier
- air permeable
- pad does not stick to the wound surface

SängerPlast Classic Cotton is the reliable textile wound plaster for rapid treatment of minor injuries. Adheres well. Easily removable from skin. Protects minor injuries from dirt. Non-sterile.



## CLASSIC COTTON

Size	Art.-No.
0.5 m x 6 cm	03 471
1 m x 6 cm	03 470

Packing: Art.-No. 03 471: 100 polybags per box  
Art.-No. 03 470: 50 polybags per box



# Plaster Cutter

Suitable for SängerPlast Pore and SängerPlast Trans, with plastic cutting edge. Flap open to place plaster roll inside. Also suitable as dispenser for the pocket. Comes without plaster rolls.



## PLASTER CUTTER

Size	Art.-No.
<b>blue</b>	
1.25 cm	02 625
2.5 cm	02 626
<b>green</b>	
1.25 cm	02 624
2.5 cm	02 627

Packing:  
1.25 cm: 500 pieces per outer carton  
2.5 cm: 250 pieces per outer carton

# Rescue Sheet

According to DIN 13 164. Protects patients from cold weather, heat and from dirt. Also avoid persons getting exposed to cold or UV-sunrays.  
Protection against cold: Silver side inside towards the body to retain body heat.  
Protection against external heat: Silver side outside to reflect heat and sun rays.



## RESCUE SHEET

Art.-No.	Size
03 620	160 x 210 cm

Packing:  
1 piece per polybag  
200 pieces per outer carton



- Non-slip, even when repeatedly re-opened or re-fitting the bandage.
- Light weight, porous and lets skin breathe.
- Soft but tension resistant material, with a tensile strength of min 85N.
- Comfortable bandage which does not require clips to deliver controlled compression.
- Protects wound bandages.
- Water repellent: Will not come off by water or sweat.

support effect:



## HAFT-O-MULL

Size	Art.-No.	Rolls/Box
4 x 400 cm	15 370	12
10 x 400 cm	15 373	12

Packing:  
Individually packed in PP-bag.  
Mixed colours per carton.  
Quantities per above table in strong, white carton with flap-lid.

## HAFT-O-DESIGN

Size	Art.-No.	Rolls/Box
5 x 450 cm	15 351	36
10 x 450 cm	15 353	18
15 x 450 cm	15 354	12

Packing:  
Individually packed in PP-bag.  
Mixed printing design and colours per carton.  
Quantities per above table in strong, white carton with flap-lid.

## HAFT-O-PREN

Size	Art.-No.	Rolls/Box
5 x 450 cm	15 341	36
10 x 450 cm	15 343	18

Packing:  
Individually packed in PP-bag.  
Mixed colours per carton.  
Quantities per above table in strong, white carton with flap-lid.



## Pre-Tape Foam Padding

Professional quality, pre-taping polyurethane foam underwrap. Protects the skin from friction and chafing, easily removable without residue, hand tearable. Breathable. Colour: white.  
Material: Polyurethane foam



## Sänger Pre-Tape Padding Bandage

Made of polyester. As underlayer underneath stiff or semi-stiff bandages, i.e. sport tape, cast or zinc-cast tape. Also suitable for padding or mechanical supporting aid of compression bandages. Soft surface structure, air permeable, skin friendly.  
Material: Polyester



## HAFT-O-FLEX

Size	Art.-No.	Rolls/Box
2.5 x 450 cm	15 320	60
5 x 450 cm	15 321	36
7.5 x 450 cm	15 322	24
10 x 450 cm	15 323	18

Packing:  
Individually packed in PP-bag.  
Mixed colours per carton.  
Quantities per above table in strong, white carton with flap-lid.

## HAFT-O-VET

Size	Art.-No.	Rolls/Box
5 x 450 cm	15 311	36
10 x 450 cm	15 313	18
15 x 450 cm	15 314	12

Packing:  
Individually packed in PP-bag.  
Mixed colours per carton.  
Quantities per above table in strong, white carton with flap-lid.

## HAFT-O-RIP

Size	Art.-No.	Rolls/Box
5 x 450 cm	15 331	36
10 x 450 cm	15 333	18

Packing:  
Individually packed in PP-bag.  
Mixed colours per carton.  
Quantities per above table in strong, white carton with flap-lid.



## PRE-TAPE

Size	Thickness	Art.-No.
7 cm x 27 m	0.8 mm	15 396

Packing:  
1 roll per Polybag  
12 rolls per box  
8 boxes per outer carton



## PRE-TAPE

Size	Thickness	Art.-No.
5 cm x 360 cm	5 - 6 mm	03 410
7.5 cm x 360 cm	5 - 6 mm	03 411
10 cm x 360 cm	5 - 6 mm	03 412

Packing: 10 pieces per Polybag  
03 410: 720 pieces per outer carton  
03 411: 480 pieces per outer carton  
03 412: 360 pieces per outer carton



## Cold-Warm Gel Pads

Cold-Warm Gel Pads for multiple use to relieve pain. Made from strong PE outer bag for long durability. Adapts well to bodyshape. Easy to clean.  
**Contents:** Non-toxic



### COLD-WARM GEL PADS

Size	Art.-No.
5.5 x 12 cm	03 631
13 x 14 cm	03 632
12 x 29 cm	03 633

**Packing:**  
 Art.-No. 03 631 100 pcs./box  
 Art.-No. 03 632 50 pcs./box  
 Art.-No. 03 633 25 pcs./box



## Ice Bag

Ideal for treating head aches and inflammation. Fill the ice bag with ice cubes or ice water and place the bag for relief. Ice bag is covered with cloth material and can be hand washed. Ice bag is flexible and fits to any body contour. Ideally suited for tooth aches or tooth extractions. For heat therapy fill the bag with hand warm water.

**Diameter:** 12 cm



### ICE BAG

Art.-No.  
03 611

**Packing:**  
1 piece per box

## Dispo Ice-Spray

Medical cooling spray for quick treatment of sports injuries. Spray from a distance of approx. 20 cm. FCKW-free.

**Content:** 200 ml / 400 ml



### DISPO ICE-SPRAY

Art.-No.	Content
03 606	200 ml
03 616	400 ml

**Packing:**  
24 Units per outer carton

Standards: Med II

## Instant Ice Bag

Instant ice wherever and whenever needed. Without using the refrigerator, deep freezer or cooling bag. Ideal for outdoor activities, in sports, for families or when travelling. Just squeeze the inner water bag and contents of instant ice bag starts its endothermic reaction obtaining an instant cooling temperature between -4°C to -8°C within seconds. Shake the bag briefly and apply to the area concerned. This will give instantaneous relief.



### INSTANT ICE BAG

Art.No  
03 630

**Packing:**  
25 pieces per outer carton

Cooling-temperature: -4°C to -8°C  
 Material of instant ice bag:  
 Polyethylene, watertight, velour surface  
 Dimensions: 14 x 18 cm

Standards: Med I

## Gel-o-Pack



Gel-o-Pack is the ideal harness to keep a cold-warm pad in the position without slipping away. Suitable for knee, wrist, neck, forehead or back.

The adjustable velcro seal closure allows flexibility. Replacement cool pads can be placed in the fridge or freezer to provide continuous cold treatment if required.  
**Contents per pack:** 1 harness with cold-warm pad  
 Replacement cold-warm pads available.

### GEL-O-PACK

Description	Art.-No.
<b>1</b> Wrist L/XL Replacement gel pads, 2 pcs. / pack	15 382 15 383
<b>2</b> Knee Replacement gel pad knee	15 392 15 393
<b>3</b> Forehead/Head/Neck/Shoulder Replacement gel pads wrist, 2 pcs. / pack 2 x Cold-Warm Pads	15 384 15 383 03 631
<b>4</b> Back Length 110 cm, extendable to 125 cm Replacement gel pads, 2 pcs. / pack	15 390 15 391

Wrist support with additional metal enforcement if required providing rigidity.





## Patient's Serviettes with Pouch

Surface layer of strong absorbing paper with under layer of liquid resistant PE-foil to protect the patient's clothing. With pouch at the bottom to additionally collect liquids. With loop around the neck for easy donning.

Size: 38 cm x 68 cm with additional 10 cm pouch.



### PATIENT'S SERVIETTES

Colour	Art.-No.
blue	03 376

**Packing:**  
600 pieces per carton



Med I

## Dental Patient Serviettes



2-layer cellulose material with water-proof underside as third layer for secure protection of patients against splashes or dirt marks on clothes. Highly absorbent. Size: 45 x 33 cm. Non-sterile.



Med I



### PATIENT SERVIETTES

Colour	Art.-No.
white	03 320
green	03 321
blue	03 322
yellow	03 323
pink	03 324
orange	03 319
Hippo	03 325

**Packing:**  
500 serviettes per pack

### DISPENSERBOX

Colour	Art.-No.
Transparent	04 321

**Measurements:** H 15 cm x W 19 cm x T 17 cm  
**Packing:**  
1 Dispenserbox per outer carton

## Patient's Serviettes Tear-Off

Surface layer of strong absorbing paper with under layer of liquid resistant PE-foil to protect the patient's clothing. For tying at the back. On roll with tear-off perforations

**Dispenser for patient's serviettes**  
white 04 322



Med I



### PATIENT'S SERVIETTES

Size: 50 x 60 cm

Colour	Art.-No.
white	03 370
blue	03 371
green	03 372

Size: 50 x 80 cm

white	03 373
blue	03 374
green	03 375

**Packing size 50 cm x 60 cm:**  
80 serviettes per roll 6 rolls per carton  
**Packing size 50 cm x 80 cm:**  
60 serviettes per roll 6 rolls per carton

## Rubber Hot Water Bottle

Made of high quality rubber with a world-wide reputation of outstanding quality and safety. 2 Litre rubber hot water bottles from Sanger's own production are CE approved. With user and cleaning instruction sheet in each carton. Many other types and sizes can be supplied. Request for our full hot water bottle catalogue. Colours: red and blue



Med I



### HOT WATER BOTTLE

Colour	Art.-No.
<b>One side ribbed</b>	
red	13 694
blue	13 693
<b>Both sides ribbed</b>	
red	13 696
blue	13 695

**Packing:**  
30 pieces per colour in carton with user instructions

## Denture Box

For hygienic storage and cleaning, for corrective dentures and replacement dentures. With strainer insert to lift the denture out of the cleaning solution and for rinsing. Closes well.

Height 55 mm x Width 95 mm x Depth 75 mm



Med I



### DENTURE BOX

Colour	Art.-No.
white	03 600
blue	03 601
turquoise	03 609

**Packing:**  
1 piece per box  
200 boxes per outer carton

## One-Way Toothbrush

with toothpaste

Ready to use, with toothpaste on the bristles. Packed individually in sealed PE sleeves. Can wet with water for easy teeth brushing. Ideally suited for dental practises, hospitals, hotels and institutions.



### TOOTHBRUSH

Art.-No.	03 608
<b>Packing:</b>	Each piece in a polybag 150 pieces per dispenser box

## Wooden Tongue Depressors

With blunted edges, has been treated for splinter free surface. Ideal for patient examinations, for application of creams and gels, for handling paste and other non-liquid substances. Non-sterile.



### TONGUE DEPRESSORS

Art.-No.	03 360
<b>Packing:</b>	100 pieces per box 100 boxes per outer carton



Med I

## Cellulose Swab

Ready to use, pre cut swabs on rolls. Size: 42 x 45 mm. Made of bleached cellulose according to DIN 19310. Absorb liquid well. Suitable as medium to apply local skin disinfectant prior to injections or infusions. Suitable to clean the skin or to be used as absorbing material. Non-sterile.



### CELLULOSE SWAB

Art.-No.	03 450
<b>Packing:</b>	2 rolls of 500 pieces per foil pack 12 pack per outer carton



Med I

## Sterilisation Foil

Rolls and bags made of polypropylene foil and paper for steam sterilisation. With sterile indicator. Materials meets the international and European standards.

### Roll sealing:

Cut to required length from roll and seal the bag on both sides with a suitable sealing machine.

### Bags:

Each bag is pre-sealed on three sides and has a self adhesive edge on the fourth side.



### STERILISATION FOIL Roll

Size	Art.-No.	Rolls/Box
50 mm x 200 m	02 445	8
75 mm x 200 m	02 446	8
100 mm x 200 m	02 447	4

### STERILISATION FOIL Bag

Size	Art.-No.	Pcs./Box
90 x 230 mm	02 444	200
300 x 450 mm	02 449	200

## Tray Paper

Ideal for standard metal trays. This filter paper can be used on the base of the tray. Weight: 55g/sqm +/-2g. Cellulose free, bleached without use of chlorine or dioxin. No recycling materials. Free of formaldehyde and optical bleaching agents. Colour fast. Available in colours white, blue, green and yellow.

Size: 28 x 18cm.



### TRAY PAPER

Colour	Art.-No.
white	03 326
blue	03 327
green	03 328
yellow	03 329

**Packing:**  
250 pieces shrink wrapped  
10 pack per outer carton

## Disposable Cups

Strong Polystyrene cups suitable for cold or warm liquids. Strong material for easy handling.



### DISPOSABLE CUPS 200 ML

Colour	Art.-No.
white	03 660

### Packing:

100 pieces per pack  
3,000 pieces per outer carton

### DISPOSABLE CUPS 175 ML

Colour	Art.-No.
white	03 667

### Packing:

100 pieces per pack  
3,000 pieces per outer carton

### DISPOSABLE CUPS 200 ML

Colour	Art.-No.
<b>200 ml</b>	
flieger	03 658
orange	03 659
light blue	03 661
pink	03 662
yellow	03 663
green	03 664
dark blue	03 665
red	03 666


### Packing:

50 pieces per pack  
1,500 pieces per outer carton

## One-Way Washing Towel

Blue and white colour non-woven material with soft surface texture. Edges sewn and extremely soft and durable. Colour: blue/white, non-sterile.

Size: 16 x 22 cm



### WASHING TOWEL

Art.-No.	Size
03 1162	16 x 22 cm

**Packing:**  
100 pieces per polybag  
10 polybags per outer carton

## One-Way Shaver

### with Twin Blades

For a quick and gentle shave of body hair of all lengths. Each shaver is closed with a protection cap. A small comb in front of the blades positions the hair in the right direction for a perfect cut. Thorough cutting action due to twin blades behind each other. The blades are made of stainless steel. Colour: blue.



### ONE-WAY SHAVER

Art.-No.
03 610

**Packing:**  
Each piece with protection cap  
100 pieces per box  
10 boxes per outer carton

## Urine Cup TAPVAL

Tightly locked container, even after numerous removals of test quantities. Save valuable working time: When removing test quantities there is no need to open or cut the lid. Just insert the test strip through the snout supplied or insert a pipette through the special TAPVAL plug. After removal, the urine cup will be as tightly closed as before. Special space for marking, 100 ml scale. Non-sterile.



### TAPVAL

Art.-No.	Volume
03 100	100 ml

**Packing:**  
600 pieces per carton



## Pill Cutter / Pill Crusher with 3 compartments



Pill cutter with blade in the lid for exact cutting of tablets, capsules, pills etc with 3 compartments for storage. The pill crusher is situated at the opposite end. Crushing is done by tightening the lid.

### PILL CUTTER / PILL CRUSHER

Size	Art.-No.
40 x 63 mm	03 614

### Packing:

1 pieces per dispenser box  
144 boxes per outer carton



Compartments



## Pill Cutter

Ideal for exact cutting of tablets, capsules and pills. Stainless blade fitted to the pill cutters lid cuts the pills in half. Two additional compartments are designed for keeping cut and uncut pills stored. Colour: Blue



### PILL CUTTER

Art.-No. 03 604

**Packing:**

1 piece per blister pack  
500 pieces per outer carton



## 7 Day Pill Box

Days are visibly marked on the cover lid. Seven separate compartments, one for each day for the required quantum of pills. Handy to carry in pocket or hand bag. When opening one day's compartment all other six compartments remain closed.



### 7 DAY PILL BOX

Art.-No. 03 602

**Packing:**

1 piece per blister bag  
400 pieces per outer carton



## Pill Organiser 7 x 4

with Braille (letters for blind people)

Extremely practical pill planer for seven days of removable pill boxes, containing four compartments for each day, stating "morning, mid-day, evening, night". Ideal pill planner for hospitals, homes for the elderly or handicapped peoples homes. Days description and Braille letters are available in German or English language. Marking: Monday till Sunday and morning, mid day, evening, night. Standard marking in German. English marking is available on request.



### PILL ORGANIZER 7X4

Art.-No. 03 603

**Packing:**

1 piece per blister pack  
100 pieces per outer carton



## Pill Crusher

For crushing of tablets, capsules, pills or other medicine. Insert the pill into the compartment and crush it by closing the lid tightly.

For save keeping collect the crushed pills into the lid mould and pour into the bottom compartment underneath milling compartment. A second compartment is available for storing uncrushed pills.

**Material:** Clear plastic with three compartments and two red colour lids.



### PILL CRUSHER

Art.-No. 03 605

**Packing:**


1 piece per box  
250 pieces per outer carton



## Disposable products for care and wellness

- One-way products
- No washing or drying
- Hygienically clean

Soft and gentle to the skin, giving a luxurious feeling, made of natural bamboo and viscose fibres for one-way use. Save cost of washing, drying and ironing. This natural material is bio-degradable. It absorbs high degree of moisture and allows the skin to breathe. Will stay in shape when wet and in use.



### WASHING GLOVE

Size: 15 x 23 cm, sewn seams

Art.-No. 04 502

**Packing:**

1 pc. / polybag  
250 pcs. per carton

### BATHING TOWEL

Size: 100 x 100 cm, sewn seams

Art.-No. 04 506 258

**Packing:**

1 pc. / polybag  
50 pcs per carton



### BATHING TOWEL

Size: 70 x 120 cm, seamless

Art.-No. 04 501

**Packing:**

1 pc. / polybag  
50 pcs per carton

### TOWEL

Size: 50 x 70 cm, seamless

Art.-No. 04 500

**Packing:**

1 pc. / polybag  
50 pcs. per carton



### FLOOR MAT

Size: 50 x 70 cm, with sewn edge


Art.-No. 04 504

**Packing:**

1 pc. / polybag  
50 pcs. per carton

## Examination Table Cover

For quick and easy covering of examination tables or chairs. No washing required, One-way product.  
**Material:** Polypropylen



### EXAM TABLE COVER WITH 2 ELASTIC LOOPS

Measurements:	Art.-No.
90 x 190 cm	04 510

**Packing:**  
130 Pieces/Carton

## Examination Table Sheets

Surface layer of strong absorbent paper with one layer of liquid resistant PE-foil underneath to protect medical examination tables from getting soiled. Keeps examination table clean enabling speedy change of sheets in between patient.

### EXAM TABLE SHEETS

Colour	Length	Art.-No.
white	50 m	03 377
blue	50 m	03 378
green	50 m	03 379
white	100 m	03 387
blue	100 m	03 388
green	100 m	03 389

**Packing:** 6 rolls per carton



## Disposable Slippers

Disposable Slippers with firm EVA sole and covered toes.  
**Size:** Length approx. 28 cm



**DISPOSABLE SLIPPERS**  
Slipper with firm EVA sole, with covered toes.

**Art.-No.** 04 508

**Packing:** 1 pair per polybag  
100 polybags per carton

## Disposable Slippers

Disposable Slippers with firm EVA sole and exposed toes.  
**Size:** approx 28 cm



**DISPOSABLE SLIPPERS**  
Slipper with firm EVA sole, with exposed toes.

**Art.-No.** 04 507

**Packing:** 1 pair per polybag  
100 polybags per carton

# Anna's Compresses

cooling - warming - caring



Awarded with Arthur Fischer Inventors  
Prize Baden-Württemberg/Germany  
2nd Place 2013



ready for use - watertight – breathable



Compresses have today regained their place as an old, traditional home remedy. Depending on their purpose, they go back to traditional water treatments or ancient herbal knowledge.

Compresses are applied as home remedy, preventive measure or care when not feeling well. They work holistically and are invigorating, soothing, cooling or stimulating - depending on additions and implementation.



Anna's Compresses are breathable and waterproof. Thus they can be applied easily and quickly in lowering body temperature by evaporative cooling, without cumbersome moisture protection in bed, which usually does not allow for the main point: evaporation.



**Contents: 1 pair**


Size 1: W 35 cm x H 24 cm  
Ref: 03 791

Size 2: W 43 cm x H 32 cm  
Ref: 03 792

Size 3: W 57 cm x H 36 cm  
Ref: 03 793



## Overcome barriers



The ideal ramp solutions for buildings, staircases or other applications where different floor height levels have to be overcome. Quick and cost effective to install, without any adhesive. For wheel chairs, wheel trolleys, carts, etc. Extremely durable and easy maintenance. Attractive finish.

- recommended for outdoor use
- recommended for indoor use

### Rubber Ramps, 3-sides sloped

Article No.	Max. height ca. mm	Width ca. mm	Depth ca. mm	Weight approx kg
■ 14 632	15	900	120	1,4
■ 14 633	25	900	200	3,3
■ 14 635	35	900	280	6,2
■ 14 636	50	900	400	12,0
■ 14 637	75	900	600	26,0
■ 14 638	100	900	800	48,0
■ 14 639	125	900	1000	67,0
■ 14 631	25	1010	200	3,5
■ 14 630	50	1140	400	12,0
■ 14 543	22	1067	203	2,6
■ 14 541	32	1067	305	7,4
■ 14 542	44	1270	457	16,0
■ 14 634	56	1524	610	28,0



3-sides sloped



wedge-type

### Rubber Ramps, wedge-type

■ 14 629	89	457	533	14,0
■ 14 628	57	914	343	9,0
■ 14 540	57	762	762	20,0
■ 14 544	4	1000	40	0,2
■ 14 545	6	1000	60	0,4
■ 14 546	8	1000	80	0,6
■ 14 547	10	1000	100	0,8
■ 14 548	12	1000	110	1,1
■ 14 549	16	1000	150	2,0

## Shower and Bath Mats



- Anti-slip surface helps prevent slipping in the bath and shower
- Effective suction cups on the underside help to keep the mat secured in place
- Perforated for good drainage
- Anti-slip properties tested
- Rinse after use for easy maintenance
- Machine washable at 40°C, using standard washing powder
- Production certified to ISO 9001:2008
- Made from natural rubber
- PVC-free

### Packing:


- Individually packed in a polybag with:
- Plastic handle
  - EAN bar-code
  - Maintenance instruction
  - English language packing material available
  - 10 pieces per carton



	Shower Mats	Bath Mats	Bath Mats
	56 x 56 cm	76 x 34 cm	96 x 34 cm
Colour	Art.-No.	Art.-No.	Art.-No.
light blue	13 719	13 739	13 759
mint green	13 722	13 742	13 762
pastel purple	13 726	13 746	13 766
white	13 720	13 740	13 760
melba	13 724	13 744	13 764
yellow	13 723	13 743	13 763
orange	13 725	13 745	13 765
red	13 714	13 734	13 754
pink	13 718	13 738	13 758
bright green	13 727	13 747	13 767
blue	13 713	13 733	13 753

### Bath Mat 76 x 34 cm & Head Rest Set

Colour	Art.-No.
light blue	13 839
mint green	13 842
pastel purple	13 846
white	13 840
melba	13 844
yellow	13 843
orange	13 845
red	13 834
pink	13 838
bright green	13 847
blue	13 833



	LATEX	VINYL	NITRILE	NEOPRENE	Sänger PRIMA® NEOTRIL
Acetic Acid (Food)	++	+	++	++	
Acetic Acid concentrated	-	-	-	+	
Acetic Acid 50%					++
Acetic Acid Anhydride	++	++	=	++	
Acetone			-	-	++
Acetonitrile					++
Acrylates		++			
Acrylic Acid	++	-	+	-	
Alcohol Based					
Desinfection Agents			++		
Alcoholic Drinks	++	++	++	++	
Aldehydes			++		
Ammonium Acetate	++		++	++	
Ammonium Carbonate	++		++	++	
Ammonium Hydroxide (30%)					++
Ammonium Nitrate				++	
Animal Fats	-	+	++	++	
Aqua Regia	-	-	+	++	
Acetic Acid, Glacial	-	-	-	+	
Battery Acid	+	+	-	++	
Benzene	-	-	-	-	
Benzyl Alcohol	-	-	++	-	
Bleaching Agents	+	+	++	++	
Borax	++	++	++	++	
Brake Oil	-	+	++	++	
Bromine			++		
Bromide	-	-	++	++	
Bromopropionic Acid	-	-	-	+	
Butoxyethanol	+	-	++	++	
Butter	-	-	++	++	
Butyl Acetate	-	-	-	+	
Butylene			++		
Calcium chloride	++	++	++	++	
Calcium Hydroxide	++	++	++	++	
Calcium Nitrate	++	++	++	++	
Calcium Phosphate	++	++	++	++	
Calcium Hypochloride	++	++	++	++	
Caustic Soda Solution 3.85 %	++				
Chlorine	-	-	-	++	
Chlorine Acetone	++	-	-	++	
Chloroform	-	-	-	-	
Chlorine Benzene	-	-	-	-	
Chromic Acid	-	-	+	-	
Citric Acid 30-70%	++	++	++	++	
Cyclohexane	-	-	-	++	
Cyclohexanol	-	-	+	++	
Diacetone Alcohol	++	-	+	++	
1,2-Dichloroethane					
Dibutyl Phthalate			++		
Dichlorethane	-	-	-	-	
Diesel Fuel	-	-	+	++	
Diethanolamine	+	++	++	++	
Diethylamine			-		
Di-Isobutyl Ketone, DIBK	-	-	-	-	+
Dimethyl Sulphoxide, DMSO	-	-	-	-	+
Dithylether					
Domestic Cleaning Agents	++	++	+	++	
Ethanolamine	-	+	++	++	
Ethanol	+	++	++	++	++
Ethidium Bromide			++		
Ethyl Acetate			-		+
Fertilizer	++	++	++	++	

	LATEX	VINYL	NITRILE	NEOPRENE	Sänger PRIMA® NEOTRIL
Formaldehyde 30%	-	-	-	++	-
Formalin Solution	+				++
Formic Acid	-	-	+	-	++
Freon 133 TF	-	-	-	-	-
Freon TMC	-	-	-	-	-
Fuel Oil	-	+	++	+	
Furfural	-	-	-	-	-
Gamma-Butyrolactone	-	-	-	-	-
Glycerin	-	-	+	++	++
Glycerophthalic Colours	-	-	-	++	-
Glycole	++	++	++	++	
Hair Bleaching Agents	++	++	++	++	
Hair Dye	++	++	++	++	
Hair Fixing Agents	++	++	++	++	
Hair Shampoo	++	++	++	++	
Hexamethyldisilazane	-	-	-	+	-
Hexane	-	-	-	++	-
Hydraulic Oil		-	+		
Hydrazine	-	++	++	++	
Hydrochloric Acid 37%			++		++
Hydrofluoric Acid (48%)					++
Hydrofluoric Acid 30%	-	-	-	-	-
Hydrogen Peroxide	++	+	++	-	
Hydrogen Peroxide	-	++	++	++	
Hydroquinone	+	+	+	+	
Isobutanol	-	-	-	+	
Isobutyl Ketone	-	-	-	-	-
Isooctane	-	-	-	++	-
Isopropyl Alcohol	-	-	-	++	+
Kerosene	-	+	++	+	
Lactic Acid 85%	++	++	++	++	
Linseed Oil	-	-	-	++	++
Lubricating Oil	-	-	-	++	-
Magnesium Oxide	++	++	++	++	
MEK					+
Methanole			-		++
Methyl		+	+		
Methylacryl Acid			-		
Methyl Acryl Ester (MAE)	-	+	+		
Methyl Ethyl Ketone			-		
Methyl Isobutyl Ketone			-		
Methylene Chloride					-
Milk And Milk Products	-	-	-	++	++
Mineral Grease	-	-	-	++	-
Monomere	-	+	+		
Naphtha	-	-	-	++	-
Naphtha VM&P	-	-	-	+	-
n-Butanole	-	-	-	-	++
Nitrobenzene			-		
n-Methyl-2-Pyrrolidone	-	-	-	-	++
n-Octanol	-	-	-	++	-
Non-Alcoholic Drinks	++	++	++	++	
Oleic Acid	+	+	++	++	
Olive Oil	-	-	-	++	++
Oxalic Acid	++	++	++	++	
Oxalic Acid, Saturated	++	++	++	++	
Palmitic Acid, saturated	-	-	-	-	++
1,2-Phenylendiamines					
Solution, saturated	++				
Paraffin Oil	-	-	-	++	-
Peanut Oil	-	-	-	++	++
Pentachlorophenol	-	-	-	++	+

	LATEX	VINYL	NITRILE	NEOPRENE	Sänger PRIMA® NEOTRIL
Pentane	-	-	-	+	-
Perchloric acid 60%	++	++	++	++	
Perchloroethylene	-	-	-	-	-
Perfume & Flavours	++	++	++	++	
Petrol/Gasoline	-	-	-	++	+
Petrolether	-	-	-	++	-
Petroleum Containing Substances	-	-	-	+	-
Phenol				++	
Phosphoric Acid	++	++	++	++	
Phosphoric Acid 85%					++
Plant Protection Agents	+	+	++	++	
Potassium Carbonate	++	++	++	++	
Potassium Hydroxide, KOH, 50%	++	++	++	++	
Potassium Nitrate	++	++	++	++	
Potassium Permanganate	++	++	++	++	
Potassium Phosphate	++	++	++	++	
Potassium Sulphate	++	++	++	++	
Poultry	-	-	-	++	++
Propylene Oxide (99%)					+
Quicklime	++	++	++	++	
Slaked Lime	++	++	++	++	
Sodium Carbonate	++		++	++	
Sodium Chloride	++		++	++	
Sodium Dicarboxylate	++	++	++	++	
Sodium Hydroxide 50%				++	++
Sodium Hydroxide, concentrated	++	+	++	++	
Sodium Hydroxide, NaOH (50%)	++	++	++	++	++
Sodium Hypochlorite	++	++	++	++	
Stoddard Solvent	-	-	-	++	-
Styrene				-	
Sulfide, Disulfide, Hyposulfide	++	++	++	++	
Sulphuric Acid (50%)					++
Tetrahydrofurfuryl Alcohol				++	
Tetrahydrofuran				-	-
Thioglycole					
Acid Solution, 20%	++				
Toluene				-	
Triethanolamine	+	+	-	+	
Turbine Oil	-	-	-	++	-
Turpentine	-	-	-	++	-
Vegetable Oil	-	-	-	++	=
Washing Powder/Detergent	++	++	++	++	++
Water Colour	++	++	++	++	++
Xylene				-	

**++** usually the glove material is very suitable for the mentioned substance  
**+** usually the glove material is suitable for the mentioned substance and will only show a mild degenerative effect  
**-** gloves made of this material have only a very limited suitability for the chemical substance  
**--** gloves made of this material are not suitable for the corresponding chemical substance  
**no Symbol** no results known  
 The above list includes information generally known in the gloves trade, for which we do not resume any responsibility or guarantees.  
 It merely serves as orientation aid. Protection against chemicals and substances is always depending on many factors, i.e. glove material, contact time, temperature etc. Disposable gloves usually provide a very limited protection against chemicals.

Account No.

\_\_\_\_\_ Contact

\_\_\_\_\_ Company Name

\_\_\_\_\_ Date

\_\_\_\_\_ Name

\_\_\_\_\_ Phone-No.

\_\_\_\_\_ Street / Number

\_\_\_\_\_ Fax-No.

\_\_\_\_\_ Country / Postal Code / Town

\_\_\_\_\_ e-mail address

\_\_\_\_\_ Delivery address, if different to above address

\_\_\_\_\_ VAT-No.

I / We order for delivery:

Description	Page	Art.-No.	Size/Colour	Quantity	Price

FO\_50\_88 10/14 Copyright © 2014 Sängers GmbH Subject to technical changes without notice.

**Sänger Fax-No.: + 49 (0) 79 35/72 24 - 199**

All details about the prices and payment terms correspond to the quotation from Sängers GmbH and their terms and conditions



